











Pull Cord Alarm

Item #	UOM
13602	1/ea

- Economical way to notify staff that a patient has moved.
- Alarm easily secures to a bed or chair with a selfcontained clip.
- Activation cord is adjustable from 18"- 36".
- When user gets up, the pin connector pulls from the alarm unit activating the signal.
- Deactivate with On/Off switch.
- Two volume settings: 97 dB and 103 dB.
- Requires one 9V battery (Included).

Tamper-Proof Magnetic Pull Cord Alarm

(Comes with Screwdriver to Remove Battery Cover)

Item #	UOM
13603	

- When magnetic pull switch is removed, alarm is activated.
- Pull switch can be activated from any direction.
- Convenient On/Off switch for easy activation and deactivation.
- Easily mounts on wheelchair or bed.
- Cord adjusts from 28"- 58".
- Alligator clip allows for attachment to patient's clothes.
- 97 103 dB.
- One 9V battery included.

specificat	ions	
Warranty		Limited 1 Year
accessori	es	
13615	Replacement Cord with Clip	1/ea
136DC	Alarm DC Adapter	1/ea

specific	ations		
Warranty	1	Limited 1 Year	
accessories			
13617	Replacement Cord, Magnet and Clip	1/ea	
13618	Replacement Plastic Alarm Caddy	1/ea	
136DC	Alarm DC Adapter	1/ea	