



## ROHO® MID PROFILE™ Single Compartment Cushion

The 3.25 in. (8.5 cm) individual air cells bring you the perfect balance immersion, protection and positioning control, redistributing pressure to decrease friction and shear. With a simple inflation adjustment, your ROHO conforms to your shape, minimizing tissue deformation and maximizing blood flow.



**Pressure injury prevention and treatment**



**Adjustable**



**Stability**



**Lightweight**



**Washable**

### Dimensions

Width: 13.25 in. to 20 in. / 33.5 cm to 50.5 cm

Depth: 13.25 in. to 20 in. / 33.5 cm to 50.5 cm

Height: 3.25 in. / 8.5 cm

Weight: Approx. 3.0 lbs. / 1.4 kg

### Weight limit

No weight limit.

However, the cushion must be properly sized to the individual.

### Limited warranty

Cushion: 24 months

Cover: 6 months



### Our evidence shows the MID PROFILE Single Compartment Cushion is ideal for people who:

- Currently have, have a history of, or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Need skin protection, immersion and envelopment
- Have normal or impaired sensation