

ROLLATORS



Rollators are a mobility aid that greatly enhance the user's autonomy. A rollator allows people of reduced strength and balance a social life – as well as exercise, improving blood circulation and general health.

The 4200LX and 4200DX are mainly for indoor use, however, their durable construction allows the occasional outdoor adventure. Made of lightweight aluminum these models are easy to fold and ideal for small spaces.

4200 4200LX WITH 6 IN DELUXE WHEELS, BLUE

4205183F 4200DX SUPER LOW, RED

4205185F 4200DX SUPER LOW, BLACK

4205223F 4200DX LOW, RED

4205225F 4200DX LOW, BLACK

The 4200WX and 4200SWX models are ideal for more active users and provide ample room between the handles while still maintaining a standard overall size.

4215181F 4200WX WIDE SUPER LOW, BLUE

4215183F 4200WX WIDE SUPER LOW, RED

4215221F 4200WX WIDE LOW, BLUE

4215223F 4200WX WIDE LOW, RED

4220223F 4200SWX SUPER WIDE LOW, RED



4200
4200LX



4205225F
4200DX
LOW



4215223F
4200WX
WIDE & LOW



4220223F
4200SWX
SUPER WIDE
LOW

For complete product
specifications & parts list,
visit
PARSONSADL.COM

Série 4200

4200SWX and 4200WX

Conçus pour les utilisateurs plus actifs, ces modèles offrent suffisamment d'espace entre les poignées tout en conservant des dimensions standards.



4200SWX



4200WX

4200DX and 4200LX

Le 4200LX et le 4200DX sont essentiellement destinés à une utilisation à l'intérieur, cependant leur conception robuste permet de les utiliser à l'extérieur occasionnellement. En aluminium léger, ces modèles sont faciles à replier et sont parfaits dans les espaces limités.



4200DX



4200LX

Accessoires



Plateau

Pratique pour poser des éléments légers.
No. de modèle : 4350



Porte-canne

Simplifie le transport d'une canne.
No. de modèle : 4345

Spécifications techniques

Description	No. de modèle	Hauteur min	Hauteur max	Hauteur du siège	Largeur	Distance entre les poignées	Diamètre des roues	Poids du produit	Capacité
4200SWX	4220223	77,5cm / 30,5po	93cm / 36,5po	55,9 cm / 22po	68,6cm / 27po	61cm / 24po	20 x 3,5cm / 7,9po x 1,4po	12,2 kg / 27lb	228 kg / 500 lb
4200WX super bas bleu	4214181	77,5cm / 30,5po	89cm / 35po	48,3cm / 19po	62,2cm / 24,5po	44,5cm / 17,5po	20 x 4cm / 7,9po x 1,6po	6,6 kg / 14,5lb	113,4 kg / 250 lb
4200WX super bas rouge	4214183								
4200WX super bas gris	4214184								
4200WX régulier, bas bleu	4214221	87,6cm / 38,5po	98cm / 38,5po	57,2cm / 22,5po	62,2cm / 24,5po	44,5cm / 17,5po	20 x 4cm / 7,9po x 1,6po	6,6 kg / 14,5lb	113,4 kg / 250 lb
4200WX régulier, bas rouge	4214223								
4200WX régulier, bas gris	4214224								
4200WX large, super bas bleu	4215181	77,5cm / 30,5	89cm / 35po	48,3cm / 19po	68,6cm / 27po	52cm / 20,5po	20 x 4cm / 7,9po x 1,6po	7kg / 15,5lbs	158,8 kg / 350 lb
4200WX large, super bas rouge	4215183								
4200WX large, super bas Gris	4215184								
4200WX large, bas bleu	4215221	87,6cm / 38,5po	98cm / 38,5po	57,2cm / 22,5po	68,6cm / 27po	52cm / 20,5po	20 x 4cm / 7,9po x 1,6po	7kg / 15,5lbs	158,8 kg / 350 lb
4200WX large, bas rouge	4215223								
4200WX large, bas gris	4215224								
4200DX super bas rouge	4205183	76,2cm / 30po	86,4cm / 34po	47,6cm / 18,75po	57,2 mm / 22,5po	42cm / 16,5po	20 x 2,5cm / 7,9po x 1po	5,4 kg / 12lb	113,4 kg / 250lb
4200DX super bas noir	4205185								
4200DX bas rouge	4205223	82,6cm / 32,5po	94cm / 37po	55,2cm / 21,75po	60cm / 23,5po	44,5cm / 17,5po	20 x 2,5cm / 7,9po x 1po	5,4 kg / 12lb	113,4 kg / 250lb
4200DX bas noir	4205225								
4200LX	4200	81,3cm / 32po	91,4cm / 36po	53,3cm / 21po	60cm / 23,5po	44,5cm / 17,5po	12,7 x 2,5cm / 5po x 1po	5kg / 11lb	113,4 kg / 250lb

Seulement distribué en Canada. Spécifications supplémentaires:

Distance entre les roues : SWX 57,2cm/22,5po, WX régulier 44,5cm / 17,5po, WX large 52cm/20,5po, DX super bas 47,6cm/18,75po, DX bas 50,2cm/19,75po, LX 43,2cm/17po

Grandeur du siège : SWX 50,8 x 35,6cm/20po x 14po, WX 43 x 31,8cm/17po x 12,5po (super bas régulier moins 2,5cm/1po), DX et LX 30,5 x 30,5cm/12po x 12po

Distance replié : SWX 30,5cm/12po, WX 34,3cm/13,5po, DX et LX 30,5cm/12po