



ROHO® HIGH PROFILE® Single Compartment Cushion

The tall, 4.25 in. (10.5 cm) individual cells envelop you in air to redistribute pressure and decrease friction and shear. With a simple inflation adjustment, your ROHO will fit your exact shape to minimize tissue deformation and maximize blood flow.



Pressure injury prevention and treatment



Positioning



Adjustable



Lightweight



Washable

Dimensions

Width: 11.5 in. to 26.75 in. / 29 cm to 68 cm

Depth: 11.5 in to 20 in. / 29 cm to 50.5 cm

Height: 4.25 in. / 10.5 cm

Weight: Approx. 3.25 lbs. / 1.5 kg

Weight limit

No weight limit. However, the cushion must be properly sized to the individual.

Limited warranty

Cushion: 24 months

Cover: 6 months



Our evidence shows the ROHO HIGH PROFILE Single Compartment Cushion is ideal for people who:

- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Require accommodation or correction of pelvic asymmetry, approx. 3 in. (7.5 cm)
- Need skin protection, maximum immersion and envelopment
- Have normal or impaired sensation
- Have a history or at risk of skin/soft tissue breakdown