

GUARDIAN® Bariatric Walkers

- Higher weight capacities keep larger patients active and moving
- Extra-wide frames give larger patients a more comfortable fit
- Two-button folding capabilities let patients receive support through narrow spaces

Compare
walkers
side-by-side
on page 35



G30754B



MDS86410XWW



G07767

Premium Bariatric 650 lb (295 kg) Capacity

Reinforced aluminum frame with steel "K" cross brace and steel cross braces on each side provide durability. Model has an extra wide frame with 22" (56 cm) between the handgrips for larger patients.

Standard Bariatric 500 lb (227 kg) Capacity

Aluminum frame with front cross brace and dual side braces provides durability. Model has an extra wide frame with 20" (51 cm) between hand grips for larger patients. Available with or without bariatric wheels.

Basic Extra Duty 400 lb (181 kg) Capacity

Lightweight aluminum frame with Guardian steel "K" cross brace. Weighs only 7.5 lbs (3.4 kg) with a 400 lb (181 kg) weight capacity. Walker has 19" (48.3 cm) between handgrips giving extra room for larger patients.

Item Number	Description	Wheels Inch (cm)	Weight Capacity Lbs (kg)	Approximate User Height Inch (cm)	Width Inside Hand Grip Inch (cm)	Height Adj. Width Inch (cm)	Folded Inch (cm)	Pkg
Premium Bariatric Folding Walkers								
G30754B	Extra Wide Bariatric	None	650 (295)	5'5"-6'4" (165-193)	22 (56)	32-39 (81-99)	4.5 (11)	2/cs
Standard Bariatric Folding Walkers								
MDS86410XW	Adult Bariatric	None	500 (227)	5'5"-6'4" (165-193)	20 (51)	31-41 (79-104)	4.5 (11)	Each
MDS86410XWW	Adult Bariatric	5 (13)	500 (227)	5'5"-6'4" (165-193)	20 (51)	33.5-43.5 (85-110)	4.5 (11)	Each
Basic Bariatric Folding Walkers								
G07767	Extra Duty, Adult	None	400 (181)	5'5"-6'4" (165-193)	19 (48)	32.25-39.25 (82-100)	4.5 (11)	Each
G07768	Extra Duty, Junior	None	400 (181)	5'-6' (152-183)	19 (48)	27.25-34.25 (69-87)	4.5 (11)	Each