



Groupe  
**Medicus**

A simple, comfortable and highly functional solution for users with mild support and positioning needs, the Invacare Absolute Wheelchair Cushion features molded, dual-firmness foam for increased comfort and support.

Gentle contour in cushion aids in centering in the seat and adds comfort  
Breathable, water-repellent mesh fabric cover dissipates heat and moisture  
No-slip bottom keeps cushion in place all day  
Molded polyurethane foam is lightweight and durable  
Shaped cushion bottom for use on taut sling upholstery

<b>Dimensions:</b>	20" W x 18" D Adductor Height: 2.7" Abductor Height: 2.3" Rear Height: 2.25"
<b>Color:</b>	Black
<b>Material:</b>	Foam
<b>Product Weight:</b>	2 lb.
<b>Warranty:</b>	Limited Warranty - 1 Year
<b>Product Weight Capacity:</b>	250 lb.