

ROHO® Toilet Seat Cushion

The adjustable, air-filled cellular neoprene surface is designed for short-term use on toilets and portable bedside commodes. It snaps easily into place and is easy to clean and disinfect.

- ✓ Quickly snaps into place
- ✓ Adjustable
- ✓ Easy to clean and disinfect
- ✓ Skin/soft tissue protection



Dimensions

Height: 2.25 in. / 5.5 cm
 Width: 7.75 in. / 20 cm
 Depth: 15.75 in. / 40 cm
 Weight: 1.25 lbs. / 0.5 kg

Weight limit

None, if properly sized
 Limited warranty
 24 months

ROHO® Shower/Commode Cushion

Using ROHO's exclusive DRY FLOATATION® Technology, the adjustable, air-filled Shower/Commode cushion provides much-needed cushioning and protection for bedside commodes and shower wheelchairs. Designed for use without a cover, the cushion is easy to clean and disinfect.

- ✓ Adjustable
- ✓ Comfortable
- ✓ Easy to clean and disinfect
- ✓ Skin/soft tissue protection



Dimensions

Height: 2.25 in. / 5.5 cm
 Width: 14.75 in. / 37.5 cm
 Depth: 14.75 in. / 37.5 cm
 Weight: 1.5 lbs. / 0.7 kg

Weight limit

None, if properly sized
 Limited warranty
 24 months