

## Toilet Safety Frame with Padded Arms

HCPCS:



FIGURE A



FIGURE B



FIGURE C



- Anodized aluminum is sturdy and lightweight
- Padded armrests provide additional comfort and support
- Comes with height-adjustable (1" increments) and width-adjustable arms (2 positions) (Figure A)
- Powder-coated aluminum bracket easily attaches frame to toilet (Figure B)
- Bracket secures safety frame in place

### PRODUCTS

### Description

RTL12000	Toilet Safety Frame,1EA,4/CS
----------	------------------------------

### SPECIFICATIONS

### Description

Adjustable Height	28" (smallest)   32" (highest)
Carton Shipping Weight	20 lbs.
Depth of Frame	18.5"
Warranty	Limited Lifetime
Weight	4 lbs.
Weight Capacity	300 lbs
Width Between Arms	19" (smallest)   23.50" (largest)