

# POWERED

## Mattress Replacement Systems



### FUSION 2K™

#### Dynamic Dispersion®

#### Mattress Replacement System

##### Manufactured by MOXI Enterprise LLC

Featuring *SelectAir*™ therapy, this system creates an optimal wound-healing environment with microclimate control.

- Convert from the Fusion 1K to the 2K in seconds without replacing the mattress. Users also can upgrade to a hospital-grade low air loss system.
- The iSense Processing™ pump includes advanced microprocessors that continually monitor a patient's immersion.
- The moisture-wicking cover keeps the patient cool and dry.
- "True" low air loss systems are used in care settings to manage microclimate of the skin and improve moisture redistribution.

**Size:** FUS2K80: 80" × 36" × 6" (203.2 cm × 91.4 cm × 15.2cm)  
FUS2K84: 84" × 36" × 6" (213.4 cm × 91.4 cm × 15.2 cm)

**System weight:** Approximately 46 lbs. (20.9 kg)

**Weight limit:** 300 lbs. (136 kg)

**Warranty:** 5-year limited warranty on the mattress; 1-year limited warranty on the cover



**U.S. Medicare Code Pending:**  
**E0277 Support Surfaces**




### SelectAir® Max

#### Low Air Loss System

##### Manufactured by MOXI Enterprise LLC

*SelectAir* systems offer same benefits as *SelectAir* plus System Lock, Auto Comfort setting and Pulse Mode which provides three levels of therapeutic alternative pressure relief.

- The firm mattress mode facilitates nursing care during dressing changes and patient repositioning.
- Upright mode provides additional support when higher elevation is required, without losing clinical effectiveness.
- The Foot Pillow cushions the foot.
- The system is divided into three sections (head, seat and feet). Integrated side air bolsters reduce the risk of entrapment.
- The zippered cover is fluid-resistant and vapor permeable.
- Provides microclimate control to keep the individual's skin dry and cool.



*SelectAir*® Standard model also available.

**Size:** SAMAXSYS: Fits most standard bed frames. Custom sizes are available.

**System weight:** Control unit is 9 lbs. (4 kg); mattress is 10 lbs. (4.5 kg)

**Weight limit:** 500 lbs. (227 kg)

**Warranty:** 5-year limited warranty on the mattress; 1-year limited warranty on the cover and control unit



**U.S. Medicare Code:**  
**E0277 Group 2 Support Surfaces**