

ADVANTAGE RAIL

The Advantage Rail is a floor mounted support rail that features our innovative pivot and lock technology for improved versatility and stability. Featuring 8 different locking positions (one at every 45 degrees), the Advantage Rail is the perfect solution for getting where you need to go in rooms that will not accommodate a floor to ceiling support pole.



SPECS - ADVANTAGE RAIL

Rail Diameter	1.5"	Contractor Recommended
Rail Height Range	30" to 38"	
Rail Length	17"	
Rail Locking Positions	8 (every 45 degrees)	
Material	Alloy Steel	
Finish	Anti-microbial powder coat white (PCW)	
Code (Weight Capacity)	AR-T (300lbs) <i>Advantage Rail</i> AR-T-HD (450lbs) <i>Advantage Rail Bariatric</i>	

SPECS - ADVANTAGE POLE

Pole Diameter	1.5"
Height	54"
Code (Weight Capacity)	AR-AP-HD (450lbs) <i>Advantage Pole</i>

PRODUCT DIMENSIONS

