

















## Deluxe Steel Drop-Arm Commode Knocked-Down Frame

Item #	UOM
11125KD-1	1/cs
11125KD-2	2/cs

- Easy-to-assemble frame.
- Easy-to-release drop-arm mechanism allows for safe lateral patient transfers to and from commode.
- Durable plastic snap-on seat and lid.
- Easy-to-clean gray powder-coated steel finish.
- Plastic armrests and back provide extra comfort and support.
- Comes with 7 qt. commode bucket with carry handle and splash shield.

specificatio	ns			
Width		19.75"	between arms   23.5" Outside	
Seat Dimensi	ions		13.75" (W) x 16.5" (D)	
Height			20" - 24"	
Weight			18 lbs.	
Carton Shipp	ing Weight (1   2)		20   39 lbs.	
Weight Limit			300 lbs.	
Warranty			Limited Lifetime	
accessories				
RTL12085	Biodegradable Sanitary Commode Liner 12/bx		Commode 12/bx, 6 bx/cs	

## Deluxe Steel Drop-Arm Commode with Padded Seat Knocked-Down Frame

Item #	UOM
11125PSKD-1	1/cs

- Easy-to-assemble without tools in 4 easy steps. (Figure A - D)
- Padded open front vinyl toilet seat.
- Easy-to-release drop-arm mechanism allows for safe lateral patient transfers to and from commode.
- Easy-to-clean gray powder-coated steel finish.
- Padded armrests and back provide extra comfort and support.
- Comes with 12 qt. commode bucket with carry handle and splash shield.

specificatio	ns	
Width		19" between arms   23.75" Outside
Seat Dimensi	ons	15.25" (W) x 16.75" (D)
Height		19" - 23"
Weight		21 lbs.
Carton Shipping Weight		23 lbs.
Weight Limit		300 lbs.
Warranty		Limited Lifetime
accessories		
RTL12085	Biodegradable Sa Liner	nitary Commode 12/bx, 6 bx/cs