



XL
HEAVY DUTY
WITH RECLINE

Aquatec *XL*

Model no. 1573862

AQUATEC XL, Heavy Duty Bath Lift - White

The Aquatec *XL* offers a weight capacity of up to 375 lbs. and features a reclining back support. The specially designed base and minimum seat height of only 2 ½" allows plenty of room for a relaxing bath.

Features

- Reinforced lifting mechanism
- Can be used with or without cover materials
- Ergonomic floating hand control
- Lithium-ion rechargeable battery pack
- Lift function will not allow the lift to descend if there is not enough battery power present to raise the user
- Can be easily folded and disassembled into two parts with one hand
- Self-releasing suction cups
- 40 degrees of recline
- Increased weight capacity

Specifications

Seat Depth	19"
Seat Width	15"
Seat Width with Side Flaps Down	28"
Overall Length	24", 33" reclined
Back Height	26"
Back Width	14"
Base Width	11"
Base Length	23"
Minimum Seat Height	2 ½"
Maximum Lift Height (Optional Risers up to 19 ¾")	17"
Weight Capacity (lbs.)	375