














# How to buy

- Measure** your client to determine the size you need. (See dimension chart for help.)
- Select** the appropriate Rifton tricycle and color.
- Choose** the accessories you need.
- Add** up the total cost.
- Order** online in minutes at Rifton.com or see the last page of this catalog for other ways to order.
- Call** us any time at **800.571.8198** if you need assistance. We're always here to help.



User dimensions (inches)	R120 small	R130 medium	R140 large
Inside leg length	17–22	22–28	25–35
 <b>Key user dimension: inside leg length</b> The extended leg should reach from seat to pedal comfortably when both feet and torso are strapped into place.  <b>Important:</b> User's weight must not exceed the maximum working load. Rifton tricycles are not intended for children with strong reflexes or poorly developed trunk balance and head control. Consult your child's therapist before ordering.			
Item dimensions (inches)	R120 small	R130 medium	R140 large
Overall height	37	40	53
Overall length	45	52	65
Overall width	27	27	31
Wheel diameter	13	15	19
Seat center to extended pedal	17–22	22–28	25–35
Seat to top of trunk support	15½–22	15½–22	15½–22
Gear ratio (approximate)	1:1	1:1	1:1
Item weight (lbs)	47	56	65
Distance between trunk laterals	8–14	8–14	8–14
Max. working load (lbs)	150	160	200

				R120 small	R130 medium	R140 large	
Required components	Rifton Tricycle frame includes self-leveling pedals with straps, backrest tube with handle and tote box for storage. Specify color: Red, blue, lime, raspberry			R112 RED <input type="checkbox"/> BLU <input type="checkbox"/> LIM <input type="checkbox"/> RAS <input type="checkbox"/>	R113 RED <input type="checkbox"/> BLU <input type="checkbox"/> LIM <input type="checkbox"/> RAS <input type="checkbox"/>	R114 RED <input type="checkbox"/> BLU <input type="checkbox"/> LIM <input type="checkbox"/> RAS <input type="checkbox"/>	
	Seat		CHOOSE	Small Large	R122 R128	R132 R138	N/A (included)
	Handlebar		CHOOSE	Conventional Loop	R126 R129	R136 R139	R146 R149
	Back support		CHOOSE	Trunk support system Backrest pad	R147 R145	R147 R145	R147 R145
	Optional accessories	Headrest (requires trunk support system)		CHOOSE	Contoured Flat	R156 R157	R156 R157
Handbrake Must be factory-installed. Cannot be added after trike is manufactured.				N/A	R125	(included)	
Communication tray Max. working load (lbs) Small: 10 Large: 15			CHOOSE	Small 8½" x 12½" Large 11½" x 13½"	K225 N/A	K225 K245	K225 K245
Abductor				R121	R131	R141	
Front guide bar				R135	R135	R135	
Rear steering bar				R123	R123	R143	
Front pulley				R144	R144	R144	
Additional backrest pad				R145	R145	R145	
Stationary stand				R155	R155	R155	