How to buy

- **1. Measure** your client to determine the size you need. (See dimension chart for help.)
- **2. Select** the appropriate Rifton tricycle and color.
- **3. Choose** the accessories you need.
- **4. Add** up the total cost.
- **5. Order** online in minutes at Rifton.com or see the last page of this catalog for other ways to order.
- **6. Call** us any time at **800.571.8198** if you need assistance. We're always here to help.



User dimensions (inches)	R120 small	R130 medium	R140 large
Inside lea lenath	17-22	22-28	25-35



and torso are strapped into place.

Important: User's weight must not exceed the maximum working load. Rifton tricycles are not intended for children with strong reflexes or poorly developed trunk balance and head control. Consult your child's therapist before ordering.

Item dimensions (inches)	R120 small	R130 medium	R140 large
Overall height	37	40	53
Overall length	45	52	65
Overall width	27	27	31
Wheel diameter	13	15	19
Seat center to extended pedal	17–22	22-28	25-35
Seat to top of trunk support	15½-22	151/2-22	151/2-22
Gear ratio (approximate)	1:1	1:1	1:1
Item weight (lbs)	47	56	65
Distance between trunk laterals	8-14	8-14	8-14
Max. working load (lbs)	150	160	200

				R120 sma	ıll	R130 me	dium	R140 larg	e
Req	Rifton Tricycle frame includes self-leveling pedals with straps, backrest tube with handle and tote box for storage. Specify color: Red, blue, lime, raspberry		R112	\$885 🗆	R113	\$1,010 🗆	R114	\$1,700 🗆	
luirec			RED□ BLU	RED□ BLU□ LIM□ RAS□		RED□ BLU□ LIM□ RAS□		RED□ BLU□ LIM□ RAS□	
Required components	Seat CH OOM	Small	R122	\$175 🗆	R132	\$175 🗆	N/A		
		Large	R128	\$310 🗆	R138	\$310 🗌	(included)		
	Handlebar C+ OSE	Conventional	R126	\$105 🗆	R136	\$115 🗌	R146	\$190	
		₹ Soc	Loop	R129	\$105 🗆	R139	\$115 🗆	R149	\$190 🗌
	Back support	Trunk support system	R147	\$410 🗆	R147	\$410	R147	\$410 🗆	
	The state of the s	SE	Backrest pad	R145	\$90 🗆	R145	\$90 □	R145	\$90 🗆
Opt	Headrest (requires trunk support system)	CHO	Contoured	R156	\$185 🗆	R156	\$185 🗆	R156	\$185 🗌
Optional accessories	(requires trunk support system)	asoose	Flat	R157	\$175 🗆	R157	\$175 🗆	R157	\$175 🗆
	Handbrake Must be factory-installed. Cannot be added after trike is manufactured.		N/A		R125	\$88	(included)		
	Communication tray Max. working load (lbs) Small: 10 Large: 15	СНО	Small 8½" x 12½"	K225	\$160 🗆	K225	\$160 □	K225	\$160 🗆
		Large 11½" x 13½"	N/A		K245	\$165 🗆	K245	\$165 🗆	
	Abductor			R121	\$77 🗆	R131	\$120 🗆	R141	\$140 🗆
	Front guide bar	**		R135	\$105 □	R135	\$105 □	R135	\$105 □
	Rear steering bar			R123	\$200 🗆	R123	\$200 🗆	R143	\$220 🗆
	Front pulley	7		R144	\$36 🗆	R144	\$36 🗆	R144	\$36 🗆
	Additional backrest pad		R145	\$90 □	R145	\$90 □	R145	\$90 □	
	Stationary stand		R155	\$95 🗆	R155	\$95 □	R155	\$95 🗆	
Tot	al			R120	\$	R130	\$	R140	\$