

## COMFORT & CARE PRODUCTS


BY VITACARE MEDICAL PRODUCTS



**PRESSURE RELIEF  
EAR CUSHION**



**VISCO "MEMORY FOAM"  
POSITIONING ROLLS**



**SEATING**



**VISCO "MEMORY FOAM" PILLOWS**



**PEDIATRIC REFLUX WEDGE POSITIONING PRODUCTS**



**EYE SURGERY RECOVERY SYSTEM**



Vitacare Medical Products is your single source for quality comfort and care medical foam products.


Our products have found wide acceptance in hospitals, institutions, clinics and homes across the country. We offer an extensive line of pressure reduction and relief mattresses, seating, pillows and positioning items as well as many specialty products.

For over 20 years, we have designed and manufactured quality products focusing on performance and user comfort.

At Vitacare Medical Products we take great pride in our quality, workmanship, service and reliability.

Visit our web-site at  
[www.vitacaremedical.com](http://www.vitacaremedical.com)  
to learn more about all of our products.

Phone: 1-800-263-9068 • Fax: 1-866-622-0986



## HEEL SOLUTIONS (PREVENTION & TREATMENT)



## VITA - DELUXE HEEL BOOT & VITA - HEEL BOOST

- Designed for patients at high risk of developing heel pressure ulcers and patients that need treatment for heel pressure ulcers
- Premium Visco "Memory Foam" wedge provides premium comfort, immersion and support while offering optimum weight redistribution for the Achilles and the calf
- Visco "Memory Foam" wedge platform is angled (down) towards the knee to reduce hyper-extension
- Firm lexan plastic base helps prevent internal and external rotation of the lower extremities
- Deluxe high density foam provides a firm outer structure that provides support and durability
- Adjustable strap closes securely with heavy duty buckle system
- Cut-outs in foam help increase air circulation
- Covered in a high quality, durable, incontinent, breathable, anti-microbial cover that is easy to clean and maintain

Note: Item# 8007 offers passive foot drop assistance



Vita - Heel Boost Item# 8006 - MSRP: \$95.00  
Important: Single patient use for stage 3 & 4



Vita - Deluxe Heel Boot Item# 8007 - MSRP: \$165.00  
Important: Single patient use for stage 3 & 4

## VITA - HEEL BOOST (UNCOVERED)


- Designed as an economical choice for the prevention of heel pressure ulcers through optimum weight redistribution
- Premium Visco "Memory Foam" wedge provides premium comfort, immersion and support while offering optimum weight redistribution for the Achilles and the calf
- Nylon base helps improve lower extremity mobility
- Deluxe high density foam provides a firm outer structure that provides support and durability
- Excellent for high volume facility use
- Designed for single patient use



Item# 8005 - MSRP: \$45.00

## ELEVATING HEEL ROLL

- Designed to prevent and treat heel pressure ulcers by offloading the heel while allowing the patient to maintain mobility (in bed)
- Perfect for patients that are in constant motion while in bed
- Rounded shape ensures consistent heel elevation during movement
- Elevates the heel at almost any angle
- Covered in a high quality, durable, incontinent, breathable, anti-microbial cover that is easy to clean and maintain
- Very easy to use and adjusts with a Velcro strap



10" outer diameter, 2 1/2" inner diameter  
Item# 8009 - MSRP: \$120.00

## ELEVATING HEEL WEDGES (SINGLE AND DOUBLE LEG DESIGNS)

- Designed for patients that do not like or cannot have anything attached to their lower extremities to offload their heels
- Elevates the heel to eliminate contact between heel and mattress to aid in the prevention and treatment of heel pressure ulcers
- Very easy to use
- Covered in a high quality, durable, incontinent, breathable, anti-microbial cover that is easy to clean and maintain


Note: Item# 8012 & Item# 8014 come with a top layer of Sunmate "Memory Foam" for premium comfort and immersion



#8011 - 36"W x 20"D x 2 1/2" - 1/2"H	MSRP: \$130.00
#8012 - 36"W x 20"D x 3 - 1"H (w/ Sunmate)	MSRP: \$180.00
#8013 - 12"W x 20"D x 2 1/2" - 1/2"H	MSRP: \$70.00
#8014 - 12"W x 20"D x 3 - 1"H (w/ Sunmate)	MSRP: \$110.00

## HALF MOON HEEL POSITIONER

- Designed for patients that do not like or cannot have anything attached to their lower extremities
- Half moon design offers a Visco "Memory Foam" upper surface and firm foam core that provides superior comfort and support while offloading the heels
- Non-slip base firmly positions patient in a fixed position
- Very easy to use
- Covered in a high quality, durable, incontinent, breathable, anti-microbial cover that is easy to clean and maintain



Item# 6275 (4"H x 8"W x 34"L) - MSRP: \$120.00

[www.vitacaremedical.com](http://www.vitacaremedical.com)

Toll Free Phone: 1-800-263-9068 / Toll Free Fax: 1-866-622-0986