

Airgo eXcursion X20 Lightweight Side-fold Rollator

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Ultra-lightweight side-fold frame designed for superior stability
- Patented ergonomic handgrips absorb vibrations, while also placing your hand at a natural angle for better brake access
- Easy to fold by just a pull of the strap, and the covered soft basket folds with the frame
- Includes a folding clip for safer lifting and transporting
- Padded seat and full contoured backrest for superior comfort
- Locks in open position for added safety
- Reflective trim for added visibility in the dark
- Easy to activate hand brakes with included parking brake function
- 8" (20 cm) comfort-grip wheels with low-profile treaded tires positioned outside the handles for increased stability
- Tool-free brake adjustment
- Reflective soft basket can be removed
- BONUS curb climber included
- AADL# C225
- ADP# MW3GM0065
- Recognized by CAOT

PRODUCTS

Description

700-920	Airgo eXcursion, Reg 20", Black, 1 EA
700-920PR	Airgo eXcursion Reg 20", Purple,1EA
700-921	Airgo eXcursion, Reg 20", Blue, 1 EA
700-922	Airgo eXcursion, Reg 20", Cranberry, 1EA

SPECIFICATIONS

Description

Handle Height	30.5" to 37"
Seat Height	20"
Seat Width	17.5"
Unit Width	25"
Warranty	Lifetime on Frame, 5 Years on Handle & Brakes, All other Parts 3 Months
Weight	15.4 lbs.
Weight Capacity	300 lbs
Weight Capacity (13046 13046-HD)	300 lbs (136 Kg)

ACCESSORIES

Description

730-942	Slow brakes sfold, Black, 1PR, 25PR/CV, PR
---------	--