

NON-POWERED
Mattress Overlays



ROHO® SOFFLEX® 2
Mattress Overlay System

A cost effective solution for different care settings while providing superior skin protection.

The SOFFLEX 2 Mattress Overlay is a non-powered, air-filled, adjustable, zoned, reactive, full support surface that is intended to conform to an individual's shape to provide skin/soft tissue and deep tissue injury protection and an environment to facilitate wound healing. The SOFFLEX 2 mattress overlay provides individualized pressure, shear and friction management through force redistribution to an individual's shape.

Size: 36" x 81.75" x 3.5" (91.5 cm x 207.5 cm x 9 cm)

Weight: 6.5 lbs. (2.9 kg)

Weight limit: 330 lbs. (150 kg)

Warranty: 18-month limited warranty



U.S. Medicare Code:
E0371 Group 2 Support Surfaces

PRODUCT FEATURES

- Adjustable to shape and weight and to changes in individual's position over time
- Promotes a comfortable environment to assist with pain management
- Accommodates other positioning devices beneath the overlay
- Segmented for use in positions other than 180°
- Hydrostatic distribution of load across the entire contact area
- Provides 3" (7.5 cm) of immersion
- Does not use electricity (no noise or heat from a pump, no issues transporting patient's bed)
- Sections are easily cleaned and can be disinfected and can be patched to extend use
- Cover is machine washable and can be disinfected
- Can be re-used with multiple individuals



ROHO® PRODIGY®
Mattress Overlay System

Provide an individualized fit in order to protect skin from breakdown.

A non-powered, air-filled, adjustable, zoned, reactive, full support surface, PRODIGY is intended to conform to an individual's shape to provide skin/soft tissue and deep tissue injury protection. Flexible, interconnected air cells can be adjusted to match the individual's body contour. Each of the three sections can be adjusted separately to address specific, changing needs.

Size: 36" x 81.75" x 3.5" (91.5 cm x 207.5 cm x 9 cm)

Weight: 9.5 lbs. (4.3 kg)

Weight Limit: 300 lbs. (136 kg)

Warranty: 12-month limited warranty



U.S. Medicare Code:
E0197 Group 1 Support Surfaces

PRODUCT FEATURES

- Adjustable to shape and weight and to changes in individual's position over time
- Promotes a comfortable environment to assist with pain management
- Zoned – Three sections each 36" x 27.5" (91.5 cm x 70 cm)
- Accommodates other positioning devices beneath the overlay
- Segmented for use in positions other than 180°
- Hydrostatic distribution of load across the entire contact area
- Provides 3" (7.5 cm) of immersion
- Does not use electricity (no noise or heat from a pump, no issues transporting patient's bed)
- Sections can be easily cleaned and disinfected and may be patched to extend use
- Cover is machine washable and can be disinfected
- Can be re-used with multiple individuals