


1. Invacare® Reduced Gap Full-length Bed Rail

Model # 6629

Durable 7/8" welded steel construction

Self-adjusting spring-loaded crossbraces snap easily into place without tools for quick installation

Compatible with 5307IVC, 531 IVC, 541 IVC




2. Invacare® Clamp-On Half Rail

Model # 6630DS

Clamp on rail, no tools required

Welded steel construction with brown vein finish and chrome on the third and fourth rail for scratch-resistant movement of the adjustment bracket

Compatible with 5307IVC, 531 IVC and 541 IVC



3. Invacare® Assist Bed Rail







Model # 6632

Heavy-duty 1" powder coated, welded steel construction

3 positions: Guard, Assist, Transfer

Moves easily out of the way with a quick release button

Compatible with 5307IVC, 531 IVC, 541 IVC and 541 IVCLOW

Model #	 Rail Length*	 Rail Height Above Deck Raised	 Rail Height Above Deck Lowered	 Width Between Rails	 Shipping Weight	 Limited Warranty
6629	55"	14.5"	4.25"	36"	27.5 lb.	1 Year
6630DS	32"	16"	3.5"	37"	24 lb.	1 Year
6632	27"	14.5"	N/A	36"	21 lb.	1 Year

*Measurement on flat bed with crossbraces installed as per instruction sheet.

Invacare beds, mattresses and rails, excluding powered support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.