


Account No. _____	Drop Ship: _____
Date: _____	Name _____
P.O. Number: _____	Address _____
Buyer: _____	City _____
Marked For: _____	State _____ Zip _____
E Mail: _____	Tel. _____ Fax _____




A (Axilla to Floor Height) _____

W (Wrist to Floor) _____

G (Groin to Floor) _____


H (Hip Width) _____

Features included in standard price




Rear drag wheel with directional stabilizer
 Leg & Foot Separator
 Height adjustable seat with spring
 Height adjustable removable back support

Height adjustable chest support
 5 wheels (1 rear & 2 front with locks)
 Height and depth adjustable handle bar

	Description	HCPCS: E8002	Size	Part No.	
	Pony, Yellow		size 0	86400	<input type="checkbox"/>
	Pony, Purple		size 1	86401	<input type="checkbox"/>
	Pony, Lilac		size 2	86402	<input type="checkbox"/>

Accessories




Hip & thigh supports, set	6"w x 3"h	size 0	86408	<input type="checkbox"/>
Hip & thigh supports, set	6"w x 4"h	size 1	86410	<input type="checkbox"/>
Hip & thigh supports, set	6"w x 4"h	size 2	86412	<input type="checkbox"/>
Minimum/Maximum Width	Size 0 6"x8¾"	Size 1 6"x8¾"	Size 2 6½"x10¾"	



Abdominal cushion	size 1	86420	<input type="checkbox"/>
Abdominal cushion	size 2	86422	<input type="checkbox"/>

Accessories



Padded Armrest
Padded Armrest
Padded Armrest

size 0 86427
size 1 86428
size 2 86429

☐
☐
☐



Vertical hand grips - variable mount

all sizes 86432

☐

Select a Head support, Swing-away Bracket, and Headrest Bracket from below



Pad Style

Head support no. 1, Black Polyurethane (cover optional) 6171550
Head support no. 2, Black Polyurethane (cover optional) 6172550
Head support no. 3, Black Polyurethane (cover optional) 6173550
Head support no. 4, Black Polyurethane (cover optional) 6174550
Anatomic head support*, Small - A = 3", B = 8" 5099001
Anatomic head support*, Medium - A = 3", B = 10¼" 5099002
Anatomic head support*, Large - A = 4¼", B = 13" 5099003

☐
☐
☐
☐
☐
☐
☐

**Includes black Ducutex cover*



Headband System (For Anatomic Head Support only) 99500
Extra headband strap (For Anatomic Head Support only) 99400-5

☐
☐



Swing-away Mounting Bracket - Height, Depth, Angle & Lateral Adjustable

Turtle bar system, short 98619-25
Turtle bar system, long 98619-26

☐
☐



Connecting rod for head support, 7.75" 1 861424-3
Connecting rod for head support, 13" 2 861424-4

☐
☐

Spare parts



Knee pad, Pony

one size 86450

☐



Push handle, Pony

one size 86416

☐





Directional stabilizer for front wheel

size 0-2 86448

☐

Measurements



	Size 0	Size 1	Size 2
	cm (inch)	cm (inch)	cm (inch)
Seat height (L)	24-38 (9¼-14¾)	38-50 (14¾-19½)	48-64 (18¾-25)
Width (G)	59 (23)	69 (27)	77 (30)
Length (H)	67 (26¼)	72 (28)	104 (40½)
Height (I)	52-69 (20¼-27)	60-82 (23½-32)	68-104 (26½-40½)
Angle range (A)	32°	32°	32°
Weight	7 kg (15.5 lb)	10 kg (22 lb)	13.5 kg (30 lb)
Max. load	30 kg (66 lb)	40 kg (88 lb)	50 kg (110 lb)

NOTES