









The Pony is a gait trainer for children with disabilities. It has a frame shaped like an arrow and a central bar that holds the chest support and the seat. The Pony is available in 3 sizes.

















Back support

support can be removed for easy access into the walker.





Height adjustable seat Allows for height adjustment to fit the user and provide optimum positioning when walking. The child is not supposed to sit on the saddle.





stabilizer for added support.






Accessories











91









Standard

As standard the Pony is delivered as featured here with back support and saddle






































	Size 0 mm	Size 1 mm	Size 2 mm
Seat height (L)	240-380	380-500	480-640
Overall width (G)	590	690	770
Length (H)	670	720	1040
Overall height (I)	520-690	600-820	680-1040
Angle range (A)	32°	32°	32°
Total mass	7 kg	10 kg	13,5 kg
Max load	30 kg	40 kg	50 kg