



R82<sup>®</sup>

# Crocodile<sup>TM</sup>





**The Crocodile is a walking aid for children and young adults to help them into a more upright position. It is placed behind the child and its lightweight design makes it easy to move. The Crocodile is available in 3 sizes.**

FRAME



270 - 690 mm



390 - 1000 mm



max. 1800 mm



max. 80 kg

### Height adjustable top frame

Easy to adjust the height without using any tools.



### Multi-adjustable handles

Positioning the grips close to the body gives better support and a better walking position.




### On/off anti-reverse stops

The anti-reverse wheel stops are used for children with less balance and control.








**Accessories are easily added or removed for the carers convenience.**






## Accessories



### River, 4 years old

Watch River walk in the Crocodile size 1. It is also available in size 2 and 3.



**1 Hand brakes**  
86834/86873

**2 Forearm support**  
86858/59/60

**3 Forearm positioner**  
869370

**4 Cross bar**  
86880

**5 Back and side support**  
8680x

**6 Hip supports**  
86819

**7 Strap for side supports**  
Sn4016

**8 Back support**  
86816

## Foldable frame

Makes storage and transporting the Crocodile easy.



**9 Back support, adjustable**  
86886

**10 Flip down seat**  
8682x, 86863

**11 Sling seat**  
8682x

**12 Handle knobs**  
86861

**13 Bag**  
86839

**14 Bag**  
Size 3, 86843

**15 Anti-tips, swing-away**  
86815 (1-2), 86815-3S (3)

**16 Attendants handle**  
86866





R82°

# Crocodile™




32°



max.50 kg



590-770 mm

670-1040  
mm

	Size 0 mm	Size 1 mm	Size 2 mm
Seat height (L)	240-380	380-500	480-640
Overall width (G)	590	690	770
Length (H)	670	720	1040
Overall height (I)	520-690	600-820	680-1040
Angle range (A)	32°	32°	32°
Total mass	7 kg	10 kg	13,5 kg
Max load	30 kg	40 kg	50 kg