



R82[®]

Mustang[™]



Mustang is a highly adaptable walking aid for children and young adults and it can be configured with anterior or posterior support to meet a wide range of walking needs. The Mustang is available in 4 sizes.

FRAME



140 - 600 cm



370 - 1250 mm



max. 1300 mm



max. 80 kg

Motivating the child to walk

The user is placed at the centre of the Mustang which provides core stability and enables the user to focus energy on walking.

Easy to access

The Mustang has been constructed with accessibility in mind. The U-shaped frame and wrap around support, allows for easy access and transfer.








Natural walking pattern

The Mustang has an easy to adjust central bar for height and angle adjustment. A slight forward leaning position, combined with built-in suspension, stimulates the natural walking pattern.

Accessories are easily added or removed for the carers convenience.



Accessories



1 Chest support
8698x0

2 Cover
8698x0-x60

3 Extension piece
869225

4 Hand bar
86936x

5 Abdominal support
869260

6 Arm rests, anterior
86935x

7 Arm prompts
8693xx

8 Forearm positioner
869370

9 Sacrum support
86925x

10 Hip support
8692xx

11 Hip supports on curved bar
869222 or 869232


12 Connecting rod
869226

Standard

You gain freedom of choice as the Mustang is delivered as standard without chest support and accessories.



13




14




15



16




17




18



19



20



13 **Seat**
8691x0

14 **Sling seat**
86957x

15 **Leg pad**
869470

16 **Guide bar**
869510

17 **Storage pouch**
869560


18 **Storage bag**
869565

19 **Leg separator**
8694x0

20 **Ankle prompts**
8694xx



Mustang[™]



| | Size 1 mm | Size 2 mm | Size 3 mm | Size 4 mm |
|--|--------------|--------------|--------------|--------------|
| Overall width (G) | 520 mm | 580 mm | 650 mm | 700 mm |
| Width between wheels (G1) | 380 mm | 450 mm | 500 mm | 560 mm |
| Length with 4" wheels (H) | 750 mm | 830 mm | - | - |
| Length with 8" wheels (H) | 820 mm | 900 mm | 1000 mm | 1100 mm |
| Height with 4" wheels (I) | 400-650 mm | 600-850 mm | - | - |
| Height with 8" wheels (I) | 490-650 mm | 600-850 mm | 700-1050 mm | 870-1300 mm |
| Seat height above floor, with 4" wheels (seat = accessory) | 170-350 mm | 240-410 mm | - | - |
| Seat height above floor, with 8" wheels (seat = accessory) | 260-440 mm | 330-500 mm | 450-680 mm | 620-840 mm |
| Width between handles (handles = accessory) | 140-350 mm | 140-350 mm | 140-600 mm | 140-600 mm |
| Center bar angle (α1) | 30° | 30° | 0-25° | 0-25° |
| Minimum turning radius | 460 mm | 550 mm | 600 mm | 650 mm |
| User height, floor to axil | 400-650 mm | 600-850 mm | 700-1050 mm | 870-1300 mm |
| Weight with 4" wheels, standard without accessories | 10 kg | 11 kg | - | - |
| Weight with 8" wheels, standard without accessories | 13 kg | 14 kg | 18 kg | 19 kg |
| Max load | 30 kg | 40 kg | 60 kg | 80 kg |
| Static stability downhill | 7° | 7° | 7° | 7° |
| Static stability uphill | 15° | 15° | 15° | 15° |
| Static stability sideways | 4,5° | 4,5° | 4,5° | 4,5° |