

Mobility - Posture Control Walkers

Kaye Posture Control Walkers are designed to promote improved alignment and increase functional walking mobility for users.

Posture Control Walkers place the support behind the user, requiring the person to adapt to a more upright posture while placing his/her center of mass within the base of support of the walker.



'B' Frame Walkers

Both 'B' Frame and 'R' Frame walkers are available with many wheel configurations and numerous optional accessories.

Walker Accessories

- Add-A-Seat
- Extensor Assist Belt
- Extensor Assist Pad
- Forearm Supports
- Guide Handles
- Leg Abductor
- Pelvic Stabilizer
- Soft Sling Support
- Suspension- Conversion Kit
- Swivel Limiters

'R' Frame = Increased Width!

'R' Frame Walkers

The 'R' Frame series is available in the three largest sizes. It provides more width to accommodate users who need additional space within the walker frame due to large orthoses or those who have a wide pelvic circumference.

Walker Accessories

- Extensor Assist Belt
- Extensor Assist Pad
- Forearm Supports
- Guide Handles
- Pelvic Stabilizer
- Soft Sling Support
- Suspension Conversion Kit
- Swivel Limiters



Model No.	Height to top of handle (floor to mid-buttocks)	Distance inside handles (depth/width)	Weight Limits (lbs/kg)
-----------	--	--	---------------------------

Posture Control 'B' Frame Walkers

W1/2B*	14.5" to 20"	9" / 13.5"	60 / 27
W1B	16.5" to 22"	9" / 13.5"	60 / 27
W2B	19" to 25"	9.5" / 13.5"	85 / 39
W3B	23" to 30.5"	10.5" / 15"	130 / 60
W4B	28.5" to 36"	13.5" / 16"	180 / 82

*Height to Top of Handle on W1/2B with rear wheels is 14½" - 18½"

Posture Control 'R' Frame Wide Walkers

R2B	19" to 25"	10.5" / 15.5"	85 / 39
R3B	23" to 30.5"	12.5" / 17.5"	130 / 60
R4B	28.5" to 36"	13.5" / 19"	180 / 82