











Universal Walker Tray with Cup Holder

Item #	UOM
10124	1/ea

- Allows personal items to be carried from room to room.
- Contains one cup holder.
- Fits most manufacturers' walkers.
- Made of easy-to-clean durable plastic.
- Easy to install.
- Size: 23" (W) x 1.5" (H) x 17" (D).

Walker Tray with Cup Holders

Item #	UOM	
10125	1/ea	

- Allows personal items to be carried from room to room.
- Contains two cup holders.
- Fits most manufacturers' 2-button walkers.
- · Easily folds down when not in use.
- Made of easy-to-clean, durable plastic.
- Easy to install.
- Size: 16" (W) x 2" (H) x 12" (D).

specifications	
Warranty	Limited 1 Year

specifications	
Warranty	Limited 1 Year