

COUCH CANE

BECAUSE YOU CAN STAND WITH CONFIDENCE





LIFT CHAIR COMPATIBLE

PRODUCT #2001













FEATURES & BENEFITS

- Ergonomic Safety Handle Makes standing and sitting easy and safe.
- Saddle Bag Pouch Four pockets keep convenient items close by.
- Adaptable Adjusts in length and height to fit either side of your favorite couch or chair. Works great with a lift chair.
- Super Grip Rubber Pads Ensure stability and protect floor from damage.

SPECIFICATIONS

Height Adjustment Handle: 34" to 40" Dimensions of Handle: $6" \times 6"$ 1.25 diameter

Base Length Adjustment: 20" to 36"

Weight: 19.5 lbs

Base Feet Measurements: 4" × 7"

Weight Limit of Handle: 250 lbs

Case of 3 Dimensions: $29"L \times 19"W \times 12"H$; 60 lbs Package Dimensions: 29" L × 19" W × 4" H; 19 lbs

WHAT IS IN THE BOX





TOOLS REQUIRED



WEIGHT CAPACITY