

LECKEY[®]

Distributed by  Mobility

Pal Classroom Seat



Posture + Attention = Learning

PAL is an attractive, robust and easy to use classroom and nursery seating system for kids of 1-12 years of age, with mild to moderate postural needs.

Simple from Start to Finish

The simplicity of **PAL** enables the therapist to assess the child's needs, instantly order their **PAL** Classroom Seat and quickly and easily set up the product when it arrives. **PAL** really is the simple, fast way to help kids join in with learning.



More Stability, Less Fatigue



Designed by the Leckey team of experienced therapists and engineers to give improved stability, which helps reduce fatigue and allows the child longer periods of concentration and fine motor activity. Available in 4 sizes and 4 colors it will help kids up to 12 years of age to participate, learn and give their full attention in the classroom.

Specifications

	Size 1	Size 2	Size 3	Size 4
Age	1-3	2-5	4-8	7-12
Seat Depth	7.1" - 9.4"	8.7" - 11.4"	9.8" - 13.0"	12" - 16"
Seat Width	5.9" - 9.8"	6.7" - 10.6"	9.1" - 13.0"	11.1" - 15.2"
Seat Weight	9.6 lbs.	13.2 lbs.	17.9 lbs.	22 lbs.
Weight Capacity	66 lbs.	77 lbs.	88 lbs.	110 lbs.
Available Colors	Pink, Green, Blue, Orange			