

Special

Application Surfaces



Size: RECLINERCUSTC: Customized to fit individual per order.
Weight: 8.5 lbs. (4.2 kg) (based on two 1R2020Cs)
Weight limit: None if properly sized to the individual
Warranty: 2-year limited warranty

ROHO® RECLINER SYSTEM

A wide range of sizes to fit most chair sizes, the DRY FLOATATION recliner system provides back and seating surface protection.

- Maintaining a person’s shape and minimizing tissue deformation is essential to the prevention and healing of pressure ulcers. DRY FLOATATION cushions allow this healing process to occur, and minimize tissue deformities.



Size: SEATTOILET: 7.75" x 15.75" x 2.25" (20 cm x 40 cm x 5.5 cm)
System weight: 1.25 lbs. (.5 kg)
Warranty: 2-year limited warranty

ROHO® TOILET SEAT CUSHION

Maximize stability and minimize potential tissue breakdowns. The toilet seat cushion is designed for standard-size toilets and portable bedside commodes.

- Multi-height cells securely position the user. The outer cells are 2.25" (5.5 cm) tall; the inner cells are 1" (2.5 cm).
- Made of impermeable, anti-static, neoprene, the toilet seat cushion is easy to clean with soap and water.




Size: THYROID: 15" x 8.5" x 5.5" (38 cm x 21.5 cm x 14 cm)
Weight: 0.5 lbs. (0.2 kg)
Weight limit: None if sized to user.
Warranty: 2-year limited warranty

ROHO® THYROID PILLOW

This pillow is designed for individuals undergoing surgical procedures or who require hyperextension of the neck. Adjustable. Assists with holding the patient’s neck in place.

- The pillow accommodates most individuals.
- This air-filled pillow eliminates disposal problems associated with foam.
- Made of impermeable, anti-static neoprene, the pillow can be disinfected between uses.



Size:
SEATCMD88LP: 14.75" x 14.75" (37.5 cm x 37.5 cm)
SEATCMD89LP: 14.75" x 16.5" (37.5 cm x 42 cm)
SEATCMD99LP: 16.5" x 18.25" (42 cm x 46.5 cm)
SEATCMD910LP: 16.5" x 16.5" (42 cm x 42 cm)
Weight: 1.5 lbs (.7 kg)
Weight limit: None if properly sized to the shower chair/commode seat.
Warranty: 2-year limited warranty

ROHO® SHOWER/ COMMODE CUSHION

Designed to fit shower chairs, the commode seat cushion offers the therapeutic benefits of a ROHO® LOW PROFILE® cushion.

- DRY FLOATATION cells minimize potential tissue breakdowns.
- Made of impermeable, anti-static neoprene, the shower/commode cushion can be disinfected between uses.