

Premium Plastic Raised, Regular/Elongated Toilet Seat, with Lock With Lock



- Designed for individuals who have difficulty sitting down or standing up from the toilet
- Wide opening in front and back for personal hygiene
- Lightweight and portable
- Standard locking mechanism ensures safety (Figure A)
- Fits most elongated toilets
- No tools required for installation
- Easy-to-clean

PRODUCTS	Description
12014	Without Arms, 1/cs

SPECIFICATIONS	Description
Seat Dimensions	17" (W) x 16.5" (D) x 4.5" (H)
Weight	4 lbs.
Carton Shipping Weight	6 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime