

R82 Crocodile - E8000

Order Form



Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

Crocodile Standard Features:

- Adjustable Rear Drag Wheel
- Lightweight Durable Aluminum Frame (sz 1 & 2)
- Stable Steel & Aluminum Frame (sz 3)
- Upholstery & Cover for Top Frame
- On/Off Anti-Reverse Wheel Stop (Prevents Unit from Rolling Backwards)
- Easy to Fold
- Solid Front & Rear Wheels
- Swivel Front Wheels with Castor Locks
- Width & Height Adjustable Rotation Handles

Size 1 and 2



Size 3



Standard Model

- ☐ Crocodile Gait Trainer, Yellow
- ☐ Crocodile Gait Trainer, Orange
- ☐ Crocodile Gait Trainer, Black
- ☐ Crocodile Gait Trainer, Metallic Grey

Size

1
2
2
3

HCPCS

E8000
E8000
E8000
E8000

Part No.

86801
86802
86852
86803

Seat Options

(A)



(B)



- ☐ Solid Seat, Yellow (A)
- ☐ Solid Seat, Orange (A)
- ☐ Solid Seat, Black (A)
- ☐ Solid Seat, Metallic Grey (A)
- ☐ Sling Seat (B)
- ☐ Sling Seat (B)

Size

1
2
2
3
1
2

HCPCS

Part No.

86820
86821
86823
86863
86824
86825

Pelvic Positioning

(A)



(B)



(C)



- ☐ Pelvic Support (7" x 4.5") (A)
- ☐ Pelvic Pad w/Lateral (3.5" x 4-3/4-9") (B)
- ☐ Pelvic Pad w/Lateral (3.5" x 7 - 11 1/2") (B)
- ☐ Hip Support, Set (C)
- ☐ Lateral Strap for Pelvic Pad

Size

1, 2
1, 2
1, 2
1, 2
1, 2

HCPCS

Part No.

86816
86808
86809
86819
905011

Trunk Positioning

(A)



(B)



(C)



- ☐ Back Support, Angle & Width Adj. (A)
- ☐ Cross Bar for Back Support (B) *
- ☐ Back Support Strap (11 1/2" x 3") (C)
- ☐ Back Support Strap (13 1/2" x 4") (C)

Size

3
3
3
3

HCPCS

Part No.

86886
86880
85490-1
85490-2

(B)* Required with size 3

		<u>Frame Options</u>	<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A)		☐ Forearm Support w/Handgrip, Set (A)	1		86858
		☐ Forearm Support w/Handgrip, Set (A)	2		86859
(B)		☐ Forearm Support w/Handgrip, Set (A)	3		86860
		☐ Padded Forearm Positioner, Each (B)*	One Size		869370
(C)		☐ Handle Knobs, Set (C)	One Size		86861
		☐ Anti-Tip Device, Swing-Away, Set (D)	1, 2		86815
(D)		☐ Anti-Tip Device, Swing-Away, Set (D)	3		86815-3S
		☐ Anti-Tip Device, Fixed, Set (E)	1, 2		86815-F
(E)		☐ Anti-Tip Device, Fixed, Set (E)	3		86815-3F
		☐ Handbrake (F) ^	1, 2		86834
(F)		☐ Handbrake (F) ^	3		86873
		☐ Attendant Handle	1, 2		86866

(B)* Need quantity 2 for a set

(F)^ When using handbrakes you will lose the anti-reverse stop

		<u>Accessories</u>	<u>Size</u>	<u>HCPCS</u>	<u>Part No.</u>
(A)		☐ Storage Bag 7" x 11" (A)	1, 2		86839
		☐ Storage Bag 11" x 14" (B)	3		86843


For shipping rates, please call 1-844-876-6245

PO Box 1739 Matthews, NC 28106 • 1-844-876-6245 • Fax: 704-882-0751

Effective: 5/1/2017



Configuration

- ☐ Crocodile
- ☐ Anti-tips
- ☐ Seat
- ☐ Hip Supports
- ☐ Arm Supports



A (Total Height) _____
 B (Shoulder Height) _____
 C (Chest Height) _____
 D (Buttock Height) _____
 E (Chest Width) _____
 F (Shoulder Width) _____
 Age _____
 Weight _____

Measurements



Crocodile

	Size 1 Inches	Size 2 Inches	Size 3 Inches
Width between handles (J)	10½ - 21½"	12½ - 23½"	13½ - 27¼"
Height - floor to handles (K)	15¼" - 27¼"	21½ - 34¼"	27½ - 39½"
Width, max (G)	24"	27"	27¼"
Length (H)	28"	32¾"	34¼"
Height - floor to top frame (I)	19 - 25¾"	25¼ - 32"	30¼"
Length, folded	28"	33½"	39½"
Height, folded	12"	10½"	11½"
Weight, standard	12 lbs.	14 lbs.	26½ lbs.
User height, maximum	44¾"	58½"	71"
User weight, maximum	66 lbs.	99 lbs.	176¼ lbs.
Turning radius	29½"	30¼"	35½"