



## AGILITY Max Contour Back System

The AGILITY Max Contour Back System delivers 7 in. (18 cm) of contour depth to provide significant posterior and lateral trunk support for optimal posture and stability. It is the only deep contour back with independent, adjustable, lateral air pads for better positioning and increased comfort.

### Options:

Foam with Air Insert

Foam with Air Insert and Lumbar Support



Significant posterior back support



Lateral trunk support and control



Stability and comfort

### Dimensions

Depth with Air Insert:

7 in. / 18 cm

Width: 13 in. to 20 in. /

33 cm to 51 cm

Height: 16 in. to 20 in. /

40.5 cm to 51 cm

Weight: 4.25 lbs. to 6.5 lbs. /

2.0 kg to 3.0 kg

### Weight limit

300 lbs. / 136 kg

Back must be properly  
sized to the individual.

### Limited warranty

Back Shell & Hardware:

5 years

Air Insert: 18 months

Cover & Foam: 18 months

Custom Air Insert: 36 months

**Our evidence shows the  
AGILITY Max Contour  
Back System is ideal for  
people who:**

- Need maximum trunk support and control
- Require significant postural support
- Have or have the potential for back or shoulder pain