




## Hälsa Bed

- » Dual posts on each side increases stability at all heights
- » Attractive, durable headboard and footboards available in Nature (Maple) and Wine (Mahogany) colours
- » Removable composite deck sections for easy cleaning
- » Weight Capacity: **375 lbs**
- » Dimensions: **36" x 80"**
- » Low Position: 10" (12¼" with large caster)
- » High Position: 25" (28" with large caster)
- » Trendelenburg/Reverse Trendelenburg
- » Rolls at any height
- » Emergency battery back-up
- » Virtually noise free with low electrical consumption DC motor
- » Cover for headboard and footboard
- » Hooks for incontinence bag
- » 3 locking casters and 1 directional caster
- » Side and foot mattress retainers
- » Receptacle for trapeze bar
- » 2 receptacles for IV holder stems
- » Disassembles into 4 sections for transport/storage
- » Ideal for home care patients
- » **10 year warranty**—Frame\*
- » **2 year warranty**—Electric\*

# permobil



Easy to  
disassemble

Large caster  
option  
available



*\*Please see Hälsa Bed manual for full warranty information.*



## Hälsa Plus Bed

- » Dual posts on each side increases stability at all heights
- » Attractive, durable headboard and footboards available in Nature (Maple) and Wine (Mahogany) colours
- » Weight Capacity: **500 lbs**
- » Width expansion: **42", 48"**
- » Dimensions: **36", 39", 42", 48" x 80"**
- » Low Position: 10" (12¼" with large caster)
- » High Position: 25" (28" with large caster)
- » Trendelenburg/Reverse Trendelenburg
- » Rolls at any height
- » Emergency battery back-up
- » Virtually noise free with low electrical consumption DC motor
- » Cover for headboard and footboard
- » Hooks for incontinence bag
- » 3 locking casters and 1 directional caster
- » Side and foot mattress retainers
- » Receptacle for trapeze bar
- » 2 receptacles for IV holder stems
- » Disassembles into 4 sections for transport/storage
- » Ideal for home care patients
- » **10 year warranty**—Frame\*
- » **2 year warranty**—Electric\*

## Why choose Hälsa?

The Hälsa Bed from Permobil is the ideal home care bed. This bed boasts incredible features such as an exceptionally low positioning option and it disassembles easily for transportation and storage.

Permobil also offers the choice of two beautiful colours for the head and foot boards - Nature (Maple) or Wine (Mahogany).

Hälsa Bed has unique qualities that separate it from other electrical beds on the market. Most home health care beds have Trendelenburg capabilities as an option, but with the Hälsa Bed, the Trendelenburg feature comes standard. Additionally, most beds don't offer the capability of moving both the head and foot sections simultaneously - the Hälsa Bed does. Interested to know more? The extensive list of features are included inside.

*With all of these remarkable features, it's clear that the Hälsa Bed is the correct choice!*

## Your Permobil contacts

### **British Columbia**

Shawn Aslani  
Tel: 604-340-2709  
Shawn.Aslani@Permobil.com

### **Alberta, Saskatchewan & Northern Territories**

Jonathan Alexander  
Tel: 780-288-1423  
Jon.Alexander@Permobil.com

### **Ontario & Manitoba**

Mike Piccin  
Tel: 519-996-6294  
Mike.Piccin@Permobil.com

### **Quebec & Atlantic Canada**

Jacques Deschenes  
Tel: 514-258-2464  
Jacques.Deschenes@Permobil.com

*We would be happy to answer any of your questions or to assist you with ordering your new Hälsa Bed.*

# Hälsa & Hälsa Plus

## Electric Adjustable Beds



Hälsa



Hälsa Plus

# permobil

Montreal | Toronto | Vancouver  
Tel: 800-265-9830 | Fax: 877-636-8944  
www.permobil.ca