

XPRESSO LITE SERIES.

The Xpresso Lite allows for a more enjoyable, independent lifestyle.

Designed for the user that requires a light-weight frame, the Xpresso Lite series features low maintenance cable brakes, a soft, center folding padded seat and featherweight tires and wheels.

FEATURES.

Innovative Design Roman arch provides support, stability & style

Seat Center folding durable waterproof covered foam cushion

Collapsible Basket Foldable & removable

Wheels 8" 40 mm wide, Solid grip / EVA

Colors Ruby Red, Azure Blue, Chrome Silver, Ebony Black, Cappuccino Brown



SPECIFICATIONS.

Mini



15 lbs
18" Seat Width
18" Seat Height
22.5" Width x 26" Length
30.5" - 33.5" Handle Height
300 lbs Weight Capacity

Regular



15 lbs
18" Seat Width
21" Seat Height
22.5" Width x 26" Length
30.5" - 35.5" Handle Height
300 lbs Weight Capacity

Tall



15 lbs
18" Seat Width
24" Seat Height
22.5" Width x 26" Length
33.5" - 40" Handle Height
350 lbs Weight Capacity

OPTIONAL ACCESSORIES.

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes
- Deluxe Back Strap Pad
- Oxygen Tank Holder



Evolution
TECHNOLOGIES™

EVOLUTION TECHNOLOGIES INC.
2530 Davies Avenue,
Port Coquitlam, BC V3C2J9

tf +1 800.556.2558
info@evolution.ca
www.evolution.ca

EVOLUTION TECHNOLOGIES USA INC.
2657 Delta Ring Road
Ferndale WA 98248

tf +1 800.839.2558
info@evolutionwalker.com
www.evolutionwalker.com