

Jazzy®

ELITE ES PORTABLE




The Jazzy® Elite ES-Portable offers a combination of great style, comfort, performance and portability and disassembles into 5 convenient pieces for your active lifestyle. The in-line, front-wheel drive design allows tight turns around corners, with large front wheels for extra absorption over obstacles.

CONVENIENT DISASSEMBLY



PERFORMANCE-MATCHED
ELECTRONICS



COMFORT HIGH-BACK SEAT
WITH HEADREST

NORTH AMERICA'S FAVORITE POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.5"
Turning Radius	25.5"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 11 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	25 lbs. each
Base Weight	48 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	16" x 16-18", 18" x 18-20", 20" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

Colors



Candy
Apple Red



Viper Blue

Features

- Easy disassembly for convenient transport
- Depth adjustable, limited recline comfort high back seat available
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Seating Options*

- Depth adjustable, limited recline comfort high back seat with solid seat pan, available sizes: 16"x16"-18", 18"x18"-20"
- Deluxe contour, limited recline high back available sizes: 18"x16"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 16"x16" or 20"x16"
- Synergy® Seat

Options

- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away leg rests
- Swing-away joystick

Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Lap belt (50", 70", 80", 90")
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger
- Flag holder

*Additional charges may apply.