



End-to-End Power Chair Solutions



	JAZZY AIR®	JAZZY® ELITE ES PORTABLE	JAZZY® SELECT 6	JAZZY® 614 HD	J6
Weight Capacity	300 lbs.	300 lbs.	300 lbs.	450 lbs.	300 lbs. <sup>14</sup>
Turning Radius <sup>1</sup>	22.5"	25.5"	21.87"	22.75"	22"
Width <sup>1</sup>	25.5" <sup>13</sup>	22.75"	23.5"	25.125"	23.19"
Length <sup>1,2</sup>	38.6"	34.75"	35.375"	35.5"	34.74"
Maximum Speed <sup>3</sup>	Up to 4 mph <sup>6</sup>	Up to 4 mph	Up to 4 mph	Up to 4.3 mph	Up to 4 mph
Ground Clearance <sup>1</sup>	2.42" at footplate	1.5"	1.75"	2.25"	2.5" (frame)
Front Wheels	6" solid casters	3" solid anti-tips	6" solid casters	6" solid casters	5" solid
Drive Wheels	10"	9" solid	10" solid	14"	10" solid (pneumatic optional)
Rear Wheels	6" solid casters	6" solid casters	6" solid casters	6" solid casters	6" solid
Drivetrain	Two-motor, in-line, mid-wheel drive	Two-motor, in-line, front-wheel drive	Two-motor, in-line, mid-wheel drive	Two-motor, in-line, Mid-Wheel 6®	Two-motor, Mid-Wheel 6®
Braking System	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical
Suspension Type	Active-Trac® (ATX)	Limited	Active-Trac® (ATX)	Active-Trac® (ATX)	Active-Trac® (ATX)
High-Back Seat (max. size)	20"W x 20"D	20"W x 20"D	20"W x 20"D	24"W x 24"D	20"W x 20"D
High-Back Seat-to-Floor Range	23" - 32.9"	21.5" - 23.5"	21.25" - 23.25"	20.5" - 23.5"	20.25" - 22.25"
Synergy® Seating	N/A	20"W x 20"D	20"W x 20"D	24"W x 24"D	20"W x 20"D
Synergy® Seat-to-Floor Range	N/A	N/A	16.5" - 18.5"	17.75" - 19.75"	16.5"-17.5"
Power Elevating Seat Option <sup>7</sup>	Standard (Power adjustable seat height)	N/A	Optional	N/A	Optional
TRU-Balance®	N/A	N/A	N/A	N/A	Yes (optional)
Specialty Controls	N/A	N/A	N/A	N/A	Yes
Standard Electronics	40A, PG VR2	40A, PG GC3	40A, PG GC3 <sup>4</sup>	50A, Dynamic Shark	70A PG VSI 6-Key 70A Q-Logic NE+ 70A Q-Logic 2 EX
Battery Size (2 required)	12 volt, U-1 or 40AH <sup>12</sup>	12 volt, U-1	12 volt, U-1	12 volt, NF-22	12 volt, U-1
Standard Battery Charger	3.5A, off-board	3A, off-board	3A, off-board	5A, off-board	5A, off-board
Per-Charge Range <sup>3,11</sup>	Up to 15.29 miles <sup>10</sup>	Up to 11 miles	Up to 15 miles	Up to 11 miles	Up to 10 miles
Battery Weight <sup>9</sup>	24.5 lbs. each	24.5 lbs. each	24.5 lbs. each	37 lbs.	24.5 lbs. each
Weight of Heaviest Piece <sup>5</sup>	100 lbs.	47.6 lbs.	100 lbs.	121.5 lbs.	100.2 lbs.
Standard Seat Weight	40.4 lbs.	43 lbs. (comfort, high-back)	43 lbs. (comfort, high-back)	51.5 lbs. (high-back)	42.5 lbs. (high-back)
Standard Seat Size	18" x 18"	18" x 18"	18" x 18"-20"	20" x 18"-20"	18" x 18"
Warranty	5-year limited on frame; 2-year limited on electronics and drive motors	5-year limited on frame; 2-year limited on electronics and drive motors	5-year limited on frame; 2-year limited on electronics and drive motors	5-year limited on frame; 2-year limited on electronics and drive motors	5-year limited on frame; 2-year limited on electronics and drive motors

- 1) Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
- 2) Without front riggings.
- 3) Varies with user weight, terrain, battery type, battery condition and tire pressure.
- 4) 40A, PG 6-key VR2 required with Power Elevating Seat Option.
- 5) Does not include seat or battery weights.
- 6) 3.5 mph while in elevated position
- 7) Not available with TRU-Balance®.
- 8) Maximum weight capacity is reduced to 400 lbs. (181 kg) on 18" deep seat configurations.
- 9) Battery weight may vary based on manufacturer.
- 10) Up to 18.4 miles with 40 Ah batteries.
- 11) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.
- 12) Optional.
- 13) Width is measured at the widest part of the seat.
- 14) 250 lbs. weight capacity with TRU-Balance® Systems.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

Beamsville, Ontario • 888-570-1113  
www.pridemobility.com

