

The **Invacare® Myon™ Active wheelchair** is a comfortable, foldable, lightweight wheelchair, that is ideally suited for everyday use. Key features of this wheelchair are increased center of gravity positioning and increased seat depth and seat width. It is a shared platform with other models in the MyOn family which means that therapists and dealers can maximize opportunities for modularity and personalized adjustments for the end user.

Highlights

- A dual crossbrace design adds frame stability, resulting in easier wheeling and better utilization of the user's body movements. Optional adaptable seat depth and seat angles to ensure that each user obtains the necessary postural stability and support.
- Almost all parts on the MyOn Active frame are interchangeable with other Myon chairs. All parts such as armrests, legrests, wheels are compatible with each other.
- Choice of an active or passive rear axle position allows the dealers and therapists to customize the centre of gravity.
- MyOn Active also comes in a Heavy Duty version with a user weight of 350 lbs, which offers flexibility in a chair to suit almost everyone.



Specifications

Overall			
Width Open	12" plus seat width		
Seat Widths	15" - 24" (20",22" and 24" are for HD frames only)		
Seat Depths	16" - 20"		
Front Seat to Floor Heights			
Rear Wheel Size	Seat Slope 0°*, Front STF		
20"	14",15",16",17"		
22"	14",15",16",17",18",19"		
24"	15",16",17",18",19",20" (3° & 5° slope only)		
*Various combinations can be done. Seat Slope 3°and 5° also available. See charts on current price list for complete details.			
Kneel to Heel Range	14" STF	15"-17" STF	18"-20" STF
80° Straight Swing In/Out	13.50"-17.00"	14.25"-19.00"	18.50" - 23.50"
90° Straight Swing In/Out	13.00"-16.50"	13.75"-18.50"	18.50" - 23.00"
Elevating Legrest Straight Swing In/Out	13.25"-16.75"	14.00"-18.75"	18.25" - 23.25"
Arm Heights			
T Arm Adj.Hght Desk/Full	8.25" - 11.25		
Flip Back Arm Fixed Hght Desk/Full	9.25"		
Flip Back Arm Adj. Hght Desk/Full	7.5" - 12.5"		
Removable Hemi Dartex	8.5" - 15.5"		
Removable Desk/Full	8.5" - 13.75"		
Back Canes	Back Cane Degree	Back Cane Heights	
Fixed Canes	0°, 7°	14", 16"/17", 20"	
Adj. Height Canes	7°	14"-18", 15"-19"	
Adj. Angle Canes	7°	14"-18", 20"	
Folding Back Canes		17",20"	
Weight	approx. 33.5 lbs		
Weight Capacity	275/350 lbs		

Invacare Canada LP
www.invacare.ca

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2010

All rights reserved.
Trademarks are identified by the symbols TM and ®.
All trademarks are owned or licensed to Invacare Corporation unless otherwise noted.

Specifications are subject to change without notification.

Invacare Canada
Form No.10-790C

