

INVACARE®
PRONTO® M91™
POWER WHEELCHAIR



The **Invacare® Pronto® M91™ Power Wheelchair** keeps life simple by offering the flexibility to meet mild to moderate positioning needs, all in a base that tackles more active daily environments. Built on the foundation of Invacare® SureStep® Suspension and true center-wheel drive technology, the Pronto M91 Wheelchair offers proven, durable construction and best-in-class technologies that provide exceptional value.

Captain's Seat Features

Semi-reclining back (standard) for additional comfort
 Choice of cloth or nylon upholstery
 Width adjustable flip-back arms (standard) with optional height adjustment
 Ergonomic seating contours provide support for longer durations of use
 Adjustable height headrest for head and neck support
 Built-in seat pouch on back
 18", 20" and 22" widths

Rehab Seat Features

Adjustable seat and back angle (ASBA) rehab seat is width and depth adjustable in 1" increments
 Angle adjustable back offers range 85° to 105°
 Round canes accept other back options for additional positioning
 Transport Tie-Down Option available*

Power Tilt Only (PTO) Plus Features

17.5" seat-to-floor height for easy access to everyday environments
 Up to 55° of tilt for pressure relief (0° - 55°, 5° - 60°)
 Standard with Invacare® MK5™ NX™ 80 controller and Invacare® SPJ™ 80 single purpose joystick
 Weight capacity of 275 lb. or optional 375 lb. with heavy-duty motor package
 Easy battery access through a simple flip back design

HCPCS Code: K0848, K0849, K0850, K0851, K0856, K0858



1. True Center-Wheel Drive
2. Invacare® SureStep® Suspension
3. Stability Lock
4. Invacare® TCD™ Traction Control Design
5. Invacare® MK5™ or MK6i™ Electronics



*TRRO includes four factory-installed brackets and wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

Invacare Corporation
 www.invacare.com

USA

One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada

570 Matheson Blvd. E.,
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

©2012

Invacare Corporation.
 All rights reserved.
 Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
 Form No. 02-083
 Rev. 09/12
 120350



Model #	Overall Dimensions	Seat Dimensions	Seat-to-Floor Height	Drive Wheels	Ground Clearance	Turning Radius	Incline
M91, M91R, M91-TS	43" L** 26" W 44.5" H	16"-24" W 16"-22" D	20"-23.5"	14"	3.25"	21.5"	9°
	Total Weight	Weight Capacity	Motor	Maximum Speed	Battery	Battery Range	Limited Warranty
	274 lb.	300/400 lb Standard 275/375 lb. Power Seating	4-Pole, 4-Pole HD	6.5 mph 4.6 mph HD	22NF	Up To 13 Miles 12.3 Miles HD	5 Years Base Frame 3 Years Seat Frame 1 Year Electronics/ Motor/Gearbox

** Footboard down

Do not operate on roads, streets or highways

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.