

# Foam Mattresses

## Transcend Series



### TRANSCEND BASIC PRESSURE REDISTRIBUTION REACTIVE MATTRESS



SP03-TBCM3576  
SP03-TBCM3580  
SP03-TBCM3584

#### Available Sizes

Transcend Basic Foam Mattress 35x76  
Transcend Basic Foam Mattress 35x80  
Transcend Basic Foam Mattress 35x84



SP03-TBCM3576PRS  
SP03-TBCM3580PRS  
SP03-TBCM3584PRS

#### Available with a 4" Firm Foam Perimeter

Transcend Basic Foam Mattress 35x76  
Transcend Basic Foam Mattress 35x80  
Transcend Basic Foam Mattress 35x84

**Warranty:** 1 yr Foam Core , 1 yr Mattress Cover

### PREVENTION UP TO STAGE II PRESSURE UCLERS

#### Product Features: Transcend Basic

Deluxe removable navy nylon cover reduces friction and shearing and is fluid proof. Bottom cover is reinforced ribbed vinyl that resists punctures and tears.



### TRANSCEND PRESSURE REDISTRIBUTION REACTIVE MATTRESS



SP03-TCM3576  
SP03-TCM3580  
SP03-TCM3584

#### Available Sizes

Transcend Foam Mattress 35x76  
Transcend Foam Mattress 35x80  
Transcend Foam Mattress 35x84



SP03-TCM3576PRS  
SP03-TCM3580PRS  
SP03-TCM3584PRS

#### Available with a 4" Firm Foam Perimeter

Transcend Foam Mattress 35x76  
Transcend Foam Mattress 35x80  
Transcend Foam Mattress 35x84

**Warranty:** 5 yr Foam Core , 1 yr Mattress Cover

### PREVENTION UP TO STAGE III PRESSURE UCLERS

#### Product Features: Transcend

Upgraded Navy Urethane Stretch cover reduces friction and shearing and is fluid proof. Bottom cover is reinforced ribbed vinyl that resists punctures and tears. 2" of Premium soft foam with our unique topper surface for enhanced comfort with pillow top foam heel section for advanced pressure redistribution. 4" of premium firm high density foam for maximum support.



### TRANSCEND ADVANCED PRESSURE REDISTRIBUTION REACTIVE



SP03-TACM3576  
SP03-TACM3580  
SP03-TACM3584

#### Available Sizes

Transcend Advanced Foam Mattress 35x76  
Transcend Advanced Foam Mattress 35x80  
Transcend Advanced Foam Mattress 35x84



SP03-TACM3576PRS  
SP03-TACM3580PRS  
SP03-TACM3584PRS

#### Available with a 4" Firm Foam Perimeter

Transcend Advanced Foam Mattress 35x76  
Transcend Advanced Foam Mattress 35x80  
Transcend Advanced Foam Mattress 35x84

**Warranty:** 7 yr Foam Core , 1 yr Mattress Cover

### PREVENTION UP TO STAGE III PRESSURE UCLERS

#### Product Features: Transcend Advanced

Upgraded Navy Urethane Stretch cover reduces friction and shearing and is fluid proof. Bottom cover is reinforced upgraded black, Non-Skid Bottom that resists punctures and tears. Qualtex™ premium foam with our unique topper surface reduces pressure and increases air flow thereby reducing skin temperature and enhancing comfort. The visco elastic foam foot and heel section further redistributes pressure. Integrated high density foam around the leading edges of the mattress help keep patients centered and provide for a more secure transfer. Medium firm high density foam center offers additional support and positioning. Premium high density firm base foam provides maximum support. The Transcend Advanced features a 7yr limited warranty.