

# Advance Bath Chair



Advance Bath Chair



Shower Trolley

One size fits all bath chairs

## Accessories



LAB/#/B Abduction Belt



LAB/1/C Head support



Bath Support

## Advance Bath Chair

The Advance Bath Chair is available in 4 sizes for infants through to young adults. The product provides the user with appropriate levels of postural support ensuring they are safe and secure at all times during bathing.

The simple and easy to use adjustments mean that the product can be set up in a variety of positions in seconds depending on the user's postural requirements.

The height adjustable legs mean that the product can be set up at the appropriate height for the parent or carer to reduce the risk of back strain or injury.



## Product Sizing

Size	1		2		3		4	
Age	1 - 5		4 - 9		8 - 14		12 - 18	
	mm	inches	mm	inches	mm	inches	mm	inches
<b>User Height</b>								
Min	750	30	950	37	1150	45	1350	53
Max	1050	41	1250	49	1450	57	1650	65
<b>Advance Bath Chair</b>	LAB/1		LAB/2		LAB/3		LAB/4	
Head support	150	6	150	6	215	8.5	215	8.5
Back support	350	14	400	16	460	18	550	22
Seat depth	225	9	275	11	335	13	400	16
Calf support	260	10	260	10	340	13	340	13
<b>Seat Height</b> 50-425mm / 2-17" <b>Max Weight</b> 72kg / 159lbs <b>Total Width</b> 430mm / 17"								
<b>Shower Trolley</b>	LSF/1							
Width	600	24						
Length:								
Min	900	36						
Max	1150	45						
<b>Seat Height from floor</b> 650-700mm / 26"- 28"								
<b>Bath Support</b>	LBS/1		LBS/2					
Head support	150	6	150	6				
Back support	600	23.6	700	27.5				
Seat depth	300	12	400	15.7				
<b>Seat Height</b> 0-150mm / 0-6" <b>Max Weight</b> 62kg / 137lbs <b>Total Width</b> 400mm / 16"								

## Colour Options

The fabric, available in seven bright colours, is rinse clean, breathable, soft to the touch and can be machine washed at 60°C.

