


INVACARE® SLINGS



DRESS TOILETING HIGH SLING



Dress Toileting High Sling

Featured belt offers secure support and positioning. Easy-to-apply Dress Toileting High Sling is suitable for seat-to-seat transfers and dressing or undressing during transfers to toilet. Dress Toileting High Sling is specially designed for users with reduced head control.

Compatible with the following lifts:

- Invacare® I-Lift® Series Lifts
- Invacare® Reliant™ 450, 600 Lifts
- Invacare® 9805 Lifts
- Invacare® Jasmine™ Lift

SIZE	A	B	C	D	Waist Belt		Weight Capacity	Model #/Material
					Min	Max		Polyester Fabric
S	34.3"	28.0"	7.9"	16.9"	24.7"	46.7"	500 lb.	2478444
M	37.4"	31.5"	9.1"	19.3"	48.4"	55.1"	500 lb.	2478445
L	40.2"	35.5"	9.1"	23.6"	52.4"	59.1"	500 lb.	2478446
XL	43.3"	42.6"	10.6"	29.9"	72.4"	80.3"	550 lb.	2484780