

Room to grow



Small chair with mini kit at its smallest dimensions

- 1 Backrest height 12½"
- 2 Trunk support width 5½"
- 3 Seat depth 7"
- 4 Seat width 7"
- 5 Armrest height above seat 5"
- 6 Seat height above footboard 7"



Small chair at its largest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 11½"
- 3 Seat depth 12"
- 4 Seat width 12"
- 5 Armrest height above seat 7½"
- 6 Seat height above footboard 12"



Small standard base
R820 Rifton Activity Chair

| User dimensions (inches) | |
|---|---|
| Height with mini kit | 32–48 28–38 |
| Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended. | |
| | |
| Item dimensions (inches) | R820 small Standard base |
| Frame width | short legs: 21 long legs: 23 |
| | short legs w/ casters: 22½ long legs w/ casters: 23½ |
| Seat height above floor | short legs: 9½–12½ long legs: 18½–21½ |
| | short legs w/ casters: 13½–16½ long legs w/ casters: 18½–21½ |
| Seat angle (tilt-in-space)–degrees | 15° forward, 15° back |
| Backrest angle–degrees | 5° forward, 20° back |
| Footboard knee angle–degrees | 45°–110° |
| Seat height above footboard | 7–12 |
| Seat width with hip guides without hip guides | 7–9 12 |
| Seat depth with mini kit | 8–12 7–11 |
| Armrest height above seat | 5–7½ |
| Trunk support width | 5½–11½ |
| Backrest height | 12½–15½ |
| Headrest height above seat | 14½–21 |
| Max. working load (lbs) | 75 |



Small hi/lo base
R830 Rifton Activity Chair

| User dimensions (inches) | |
|---|-----------------------|
| Height with mini kit | 32–48 28–38 |
| Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended. | |
| | |
| Item dimensions (inches) | R830 small Hi/lo base |
| Frame width | 26 |
| Seat height above floor | 10–25 |
| Seat angle (tilt-in-space)–degrees | 15° forward, 25° back |
| Backrest angle–degrees | 5° forward, 20° back |
| Footboard knee angle–degrees | 45°–110° |
| Seat height above footboard | 7–12 |
| Seat width with hip guides without hip guides | 7–9 12 |
| Seat depth with mini kit | 8–12 7–11 |
| Armrest height above seat | 5–7½ |
| Trunk support width | 5½–11½ |
| Backrest height | 12½–15½ |
| Headrest height above seat | 14½–21 |
| Max. working load (lbs) | 75 |

Retail Order Form - January 2020

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.



| | | Standard R820 small | | Hi/lo R830 small |
|--|--------|---------------------------|--------------------------|-------------------------------|
| * Seat and back includes seatbelt (Hi/lo includes footboard and push handles) | CHOOSE | Adjustable | <input type="checkbox"/> | R821 <input type="checkbox"/> |
| | | Adjustable w/ back spring | <input type="checkbox"/> | R831 <input type="checkbox"/> |
| * Base | CHOOSE | Without spring | <input type="checkbox"/> | R822 <input type="checkbox"/> |
| | | With spring | <input type="checkbox"/> | R832 <input type="checkbox"/> |
| * Arm supports (pair) | CHOOSE | Armrests | <input type="checkbox"/> | R823 <input type="checkbox"/> |
| | | Forearm prompts | <input type="checkbox"/> | R833 <input type="checkbox"/> |
| * Adjustable legs | CHOOSE | Short | <input type="checkbox"/> | R824 <input type="checkbox"/> |
| | | Long | <input type="checkbox"/> | N/A |
| | | Short with casters | <input type="checkbox"/> | N/A |
| | | Long with casters | <input type="checkbox"/> | N/A |
| * Pads | CHOOSE | Tan | <input type="checkbox"/> | R827 <input type="checkbox"/> |
| | | Blue | <input type="checkbox"/> | R809 <input type="checkbox"/> |
| | | Red | <input type="checkbox"/> | R827 <input type="checkbox"/> |
| | | Pink | <input type="checkbox"/> | R809 <input type="checkbox"/> |
| | | Purple | <input type="checkbox"/> | R827 <input type="checkbox"/> |
| | | Green | <input type="checkbox"/> | R809 <input type="checkbox"/> |
| Headrest | CHOOSE | Contoured | <input type="checkbox"/> | R827 <input type="checkbox"/> |
| | | Flat | <input type="checkbox"/> | R809 <input type="checkbox"/> |
| | | Adjustable winged | <input type="checkbox"/> | R827 <input type="checkbox"/> |
| Laterals (pair) | CHOOSE | Small | <input type="checkbox"/> | R827 <input type="checkbox"/> |
| | | Small with strap | <input type="checkbox"/> | R809 <input type="checkbox"/> |
| Push handles (recommended for all chairs with wheels) | | | <input type="checkbox"/> | Included |
| Tray (requires armrests) | | | <input type="checkbox"/> | R667 <input type="checkbox"/> |
| Handhold | | | <input type="checkbox"/> | R649 <input type="checkbox"/> |
| Chest strap (wide) | | | <input type="checkbox"/> | R802 <input type="checkbox"/> |
| Butterfly harness | | | <input type="checkbox"/> | R834 <input type="checkbox"/> |
| Thigh belt | | | <input type="checkbox"/> | R811 <input type="checkbox"/> |
| Pelvic harness | | | <input type="checkbox"/> | R845 <input type="checkbox"/> |
| Hip guides (pair) | | | <input type="checkbox"/> | R814 <input type="checkbox"/> |
| Abductor | | | <input type="checkbox"/> | R884 <input type="checkbox"/> |
| Adductors (pair) | | | <input type="checkbox"/> | R892 <input type="checkbox"/> |
| Leg prompt | | | <input type="checkbox"/> | R854 <input type="checkbox"/> |
| Footboard (recommended for chairs with long legs) | | | <input type="checkbox"/> | Included |
| Mini kit (backrest insert and trunk support) | | | <input type="checkbox"/> | R812 <input type="checkbox"/> |
| Ankle straps (pair) | | | <input type="checkbox"/> | R899 <input type="checkbox"/> |
| Sandals (pair) | | | <input type="checkbox"/> | R944 <input type="checkbox"/> |
| Wedges (pair, requires sandals) | | | <input type="checkbox"/> | R697 <input type="checkbox"/> |
| Lumbar and seat support kit | | | <input type="checkbox"/> | R886 <input type="checkbox"/> |
| Backrest filler pad | | | <input type="checkbox"/> | R849 <input type="checkbox"/> |
| Whitmyer headrest adapter | | | <input type="checkbox"/> | R804 <input type="checkbox"/> |
| | | R820 | | R830 |

* Required

Room to grow



Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"



Medium chair at its largest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"




Medium standard base
R840 Rifton Activity Chair



Medium hi/lo base
R850 Rifton Activity Chair

User dimensions (inches)

Height 40–56



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.


Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.

Item dimensions (inches)

| R840 medium Standard base | |
|--|---|
| Frame width | short legs: 23¾ long legs: 24½ |
| | short legs w/ casters: 24½ long legs w/ casters: 25½ |
| Seat height above floor | short legs: 12½–16½ long legs: 19–23 |
| | short legs w/ casters: 16–20 long legs w/ casters: 19–23 |
| Seat angle (tilt-in-space)–degrees | 15° forward, 15° back |
| Backrest angle–degrees | 5° forward, 20° back |
| Footboard knee angle–degrees | 45°–110° |
| Seat height above footboard | 7¾–15¼ |
| Seat width with hip guides without hip guides | 8½–11½ 14 |
| Seat depth | 11–16 |
| Armrest height above seat | 7–10½ |
| Distance between laterals | 6½–12 |
| Backrest height | 15½–19 |
| Headrest height above seat | 17½–24½ |
| Max. working load (lbs) | 150 |

User dimensions (inches)

Height 40–56



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.

Item dimensions (inches)

| R850 medium Hi/lo base | |
|--|-----------------------|
| Frame width | 27½ |
| Seat height above floor | 12–23 |
| Seat angle (tilt-in-space)–degrees | 15° forward, 25° back |
| Backrest angle–degrees | 5° forward, 20° back |
| Footboard knee angle–degrees | 45°–110° |
| Seat height above footboard | 7¾–15¼ |
| Seat width with hip guides without hip guides | 8½–11½ 14 |
| Seat depth | 11–16 |
| Armrest height above seat | 7–10½ |
| Distance between laterals | 6½–12 |
| Backrest height | 15½–19 |
| Headrest height above seat | 17½–24½ |
| Max. working load (lbs) | 150 |

Retail Order Form - January 2020

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.



| | | Standard R840 medium | Hi/lo R850 medium |
|--|---------------------------|-------------------------------|-------------------------------|
| * Seat and back includes seatbelt (Hi/lo includes footboard and push handles) | CHOOSE Adjustable | R841 <input type="checkbox"/> | R851 <input type="checkbox"/> |
| | Adjustable w/ back spring | R842 <input type="checkbox"/> | R852 <input type="checkbox"/> |
| * Base | CHOOSE Without spring | R843 <input type="checkbox"/> | R853 <input type="checkbox"/> |
| | With spring | R844 <input type="checkbox"/> | N/A |
| * Arm supports (pair) | CHOOSE Armrests | R868 <input type="checkbox"/> | R868 <input type="checkbox"/> |
| | Forearm prompts | R869 <input type="checkbox"/> | R869 <input type="checkbox"/> |
| * Adjustable legs | CHOOSE Short | R893 <input type="checkbox"/> | N/A |
| | Long | R894 <input type="checkbox"/> | N/A |
| | Short with casters | R895 <input type="checkbox"/> | N/A |
| | Long with casters | R896 <input type="checkbox"/> | N/A |
| * Pads | CHOOSE Tan | R835 <input type="checkbox"/> | R835 <input type="checkbox"/> |
| | Blue | R836 <input type="checkbox"/> | R836 <input type="checkbox"/> |
| | Red | R837 <input type="checkbox"/> | R837 <input type="checkbox"/> |
| | Pink | R974 <input type="checkbox"/> | R974 <input type="checkbox"/> |
| | Purple | R975 <input type="checkbox"/> | R975 <input type="checkbox"/> |
| | Green | R976 <input type="checkbox"/> | R976 <input type="checkbox"/> |
| Headrest | CHOOSE Contoured | R897 <input type="checkbox"/> | R897 <input type="checkbox"/> |
| | Flat | R898 <input type="checkbox"/> | R898 <input type="checkbox"/> |
| | Adjustable winged | R881 <input type="checkbox"/> | R881 <input type="checkbox"/> |
| Laterals (pair) | CHOOSE Small | R828 <input type="checkbox"/> | R828 <input type="checkbox"/> |
| | Small with strap | R829 <input type="checkbox"/> | R829 <input type="checkbox"/> |
| | Large | R838 <input type="checkbox"/> | R838 <input type="checkbox"/> |
| Push handles (recommended for all chairs with wheels) | | R878 <input type="checkbox"/> | Included |
| Tray (requires armrests) | | R668 <input type="checkbox"/> | R668 <input type="checkbox"/> |
| Handhold | | R649 <input type="checkbox"/> | R649 <input type="checkbox"/> |
| Chest strap (wide) | | R863 <input type="checkbox"/> | R863 <input type="checkbox"/> |
| Butterfly harness | CHOOSE Regular | R864 <input type="checkbox"/> | R864 <input type="checkbox"/> |
| | Slim cut | R882 <input type="checkbox"/> | R882 <input type="checkbox"/> |
| Thigh belt | | R865 <input type="checkbox"/> | R865 <input type="checkbox"/> |
| Pelvic harness | | R866 <input type="checkbox"/> | R866 <input type="checkbox"/> |
| Hip guides (pair) | | R867 <input type="checkbox"/> | R867 <input type="checkbox"/> |
| Abductor | CHOOSE Small | R884 <input type="checkbox"/> | R884 <input type="checkbox"/> |
| | Large | R891 <input type="checkbox"/> | R891 <input type="checkbox"/> |
| Adductors (pair) | | R892 <input type="checkbox"/> | R892 <input type="checkbox"/> |
| Leg prompt | | R858 <input type="checkbox"/> | R858 <input type="checkbox"/> |
| Footboard (recommended for chairs with long legs) | | R848 <input type="checkbox"/> | Included |
| Ankle straps (pair) | | R899 <input type="checkbox"/> | R899 <input type="checkbox"/> |
| Sandals (pair) | | R954 <input type="checkbox"/> | R954 <input type="checkbox"/> |
| Wedges (pair, requires sandals) | | R698 <input type="checkbox"/> | R698 <input type="checkbox"/> |
| Lumbar and seat support kit | | R887 <input type="checkbox"/> | R887 <input type="checkbox"/> |
| Backrest filler pad | | R849 <input type="checkbox"/> | R849 <input type="checkbox"/> |
| Whitmyer headrest adapter | | R804 <input type="checkbox"/> | R804 <input type="checkbox"/> |
| | | R840 | R850 |

* Required

Room to grow



Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 11¼"



Large chair at its largest dimensions

- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 21¼"



Large standard base
R860 Rifton Activity Chair



Large hi/lo base
R870 Rifton Activity Chair

User dimensions (inches)

Height 50–74



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.

Item dimensions (inches)

| R860 large Standard base | |
|--|---|
| Frame width | short legs: 26¾ long legs: 26¾ |
| | short legs w/ casters: 27¾ long legs w/ casters: 27¾ |
| Seat height above floor | short legs: 16½–20½ long legs: 20–24 |
| | short legs w/ casters: 16½–20½ long legs w/ casters: 20–24 |
| Seat angle (tilt-in-space)—degrees | 15° forward, 15° back |
| Backrest angle—degrees | 5° forward, 20° back |
| Footboard knee angle—degrees | 45°–110° |
| Seat height above footboard | 11¼–21¼ |
| Seat width with hip guides without hip guides | 11–14 17 |
| Seat depth | 15–20 |
| Armrest height above seat | 7–10½ |
| Distance between laterals | 9½–14 |
| Backrest height | 19–24 |
| Headrest height above seat | 19½–29½ |
| Max. working load (lbs) | 250 |

User dimensions (inches)

Height 50–74



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.

Item dimensions (inches)

| R870 large Hi/lo base | |
|--|-----------------------|
| Frame width | 29¾ |
| Seat height above floor | 16–24 |
| Seat angle (tilt-in-space)—degrees | 15° forward, 25° back |
| Backrest angle—degrees | 5° forward, 20° back |
| Footboard knee angle—degrees | 45°–110° |
| Seat height above footboard | 11¼–21¼ |
| Seat width with hip guides without hip guides | 11–14 17 |
| Seat depth | 15–20 |
| Armrest height above seat | 7–10½ |
| Distance between laterals | 9½–14 |
| Backrest height | 19–24 |
| Headrest height above seat | 19½–29½ |
| Max. working load (lbs) | 250 |

Retail Order Form - January 2020

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.



| | | Standard R860 large | Hi/lo R870 large |
|--|---------------------------|-------------------------------|-------------------------------|
| * Seat and back includes seatbelt (Hi/lo includes footboard and push handles) | CHOOSE Adjustable | R861 <input type="checkbox"/> | R871 <input type="checkbox"/> |
| | Adjustable w/ back spring | R862 <input type="checkbox"/> | R872 <input type="checkbox"/> |
| * Base | CHOOSE Without spring | R873 <input type="checkbox"/> | R883 <input type="checkbox"/> |
| | With spring | R874 <input type="checkbox"/> | N/A |
| * Arm supports (pair) | CHOOSE Armrests | R868 <input type="checkbox"/> | R868 <input type="checkbox"/> |
| | Forearm prompts | R879 <input type="checkbox"/> | R879 <input type="checkbox"/> |
| * Adjustable legs | CHOOSE Short | R808 <input type="checkbox"/> | N/A |
| | Long | R894 <input type="checkbox"/> | N/A |
| | Short with casters | R895 <input type="checkbox"/> | N/A |
| | Long with casters | R896 <input type="checkbox"/> | N/A |
| * Pads | CHOOSE Tan | R875 <input type="checkbox"/> | R875 <input type="checkbox"/> |
| | Blue | R876 <input type="checkbox"/> | R876 <input type="checkbox"/> |
| | Red | R877 <input type="checkbox"/> | R877 <input type="checkbox"/> |
| | Pink | R977 <input type="checkbox"/> | R977 <input type="checkbox"/> |
| | Purple | R978 <input type="checkbox"/> | R978 <input type="checkbox"/> |
| | Green | R979 <input type="checkbox"/> | R979 <input type="checkbox"/> |
| Headrest | CHOOSE Contoured | R897 <input type="checkbox"/> | R897 <input type="checkbox"/> |
| | Flat | R898 <input type="checkbox"/> | R898 <input type="checkbox"/> |
| | Adjustable winged | R881 <input type="checkbox"/> | R881 <input type="checkbox"/> |
| Laterals (pair) | CHOOSE Small | R828 <input type="checkbox"/> | R828 <input type="checkbox"/> |
| | Large | R838 <input type="checkbox"/> | R838 <input type="checkbox"/> |
| | Large with strap | R839 <input type="checkbox"/> | R839 <input type="checkbox"/> |
| Push handles (recommended for all chairs with wheels) | | R878 <input type="checkbox"/> | Included |
| Tray (requires armrests) | | R669 <input type="checkbox"/> | R669 <input type="checkbox"/> |
| Handhold | | R649 <input type="checkbox"/> | R649 <input type="checkbox"/> |
| Chest strap (wide) | | R807 <input type="checkbox"/> | R807 <input type="checkbox"/> |
| Butterfly harness | CHOOSE Regular | R889 <input type="checkbox"/> | R889 <input type="checkbox"/> |
| | Slim cut | R826 <input type="checkbox"/> | R826 <input type="checkbox"/> |
| Thigh belt | | R813 <input type="checkbox"/> | R813 <input type="checkbox"/> |
| Pelvic harness | | R846 <input type="checkbox"/> | R846 <input type="checkbox"/> |
| Hip guides (pair) | | R867 <input type="checkbox"/> | R867 <input type="checkbox"/> |
| Abductor | CHOOSE Small | R884 <input type="checkbox"/> | R884 <input type="checkbox"/> |
| | Large | R891 <input type="checkbox"/> | R891 <input type="checkbox"/> |
| Adductors (pair) | | R892 <input type="checkbox"/> | R892 <input type="checkbox"/> |
| Leg prompt | | R859 <input type="checkbox"/> | R859 <input type="checkbox"/> |
| Footboard (recommended for chairs with long legs) | | R857 <input type="checkbox"/> | Included |
| Ankle straps (pair) | | R855 <input type="checkbox"/> | R855 <input type="checkbox"/> |
| Sandals (pair) | | R964 <input type="checkbox"/> | R964 <input type="checkbox"/> |
| Wedges (pair, requires sandals) | | R699 <input type="checkbox"/> | R699 <input type="checkbox"/> |
| Lumbar and seat support kit | | R888 <input type="checkbox"/> | R888 <input type="checkbox"/> |
| Backrest filler pad | | R856 <input type="checkbox"/> | R856 <input type="checkbox"/> |
| Whitmyer headrest adapter | | R804 <input type="checkbox"/> | R804 <input type="checkbox"/> |
| | | R860 | R870 |

* Required