






it is positioned above the knees, it hinders the adduction of the lower limbs. It guides the step pattern making the right sequence of movements easier while walking. It is adjustable in height and circumference to fit in the different body builds. Furthermore the extent of the stride can be adjusted. Thanks to its innovative design, the thigh straps offer three degrees of freedom: lengthwise, vertically and rotationally. They go along with the natural physiological movements of walking.




Grillo Mini Accessories



Measurements in Inches

Grillo	Mini	Small	Medium	Large
Max Capacity	77 lbs	99 lbs	198 lbs	242 lbs
User Height	29 - 35	35 - 53	43 - 65	59 - 75
Unit Weight	28 lbs	36 lbs	43 lbs	51 lbs
A Support Circumference	21 - 31	23.5 - 34	29 - 47	33 - 55
B Chest Support	21 - 30.5	27 - 39	30 - 49	32 - 56
C Pelvic Support	16.5 - 24	19 - 33	25 - 46	25 - 46
D Overall Length	28.7	29.5	36	43.5
E Overall Width	22.5	25.5	30	31.5





7 versions
 2 models: front drive and posterior,
 just one frame
folding frame
 gait trainer with constant center of gravity
 precise and continuous adjustments
 interchangeable components
 3 sizes that cover 5 sizes
 2 ergonomic and independent
 supports for trunk and pelvis

Versatile and Adjustable wide range of interchangeable components




Grillo is a gait trainer designed to help and/or make the child's and adult's self-movement easier.

Its original design and its patented height adjustment system make it unique in safety and elegance.



Its wide range of adjustments and supports make it easy-fitting to the different needs and features of the User.



Folding seat height adjustable. It can be used on posterior Grillo sizes OneTwo and ThreeFour, both on pelvic support and lumbar support. it offers a secure support in case of need, it can be folded while walking



Folding saddle height adjustable. It can be used with front drive Grillo on pelvic support and with posterior Grillo on lumbar support



(only for Grillo size Three-Four and Five)
 Pair of fixed Ø200 wheels, with pedals for parking brakes + forks.
 The frame width is reduced :
 - size ThreeFour going from cm. 76 to cm. 68,5
 - size Five going from cm. 80 to cm. 72.5

