

# R82 Manatee

## Order Form



Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_

Date: \_\_\_\_\_

Contact: \_\_\_\_\_

Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### Manatee Standard Features:

- Size 0, 1, 2 or 3
- Head/Lateral Support
- Height Adjustable Foldable Frame
- Cover in Red or Blue Soft Nylon
- Belt for Waist/Hip
- HCPCS n/a for Manatee



### Standard Model

- ☐ Manatee, Blue
- ☐ Manatee, Blue
- ☐ Manatee, Blue
- ☐ Manatee, Blue
- ☐ Manatee, Red
- ☐ Manatee, Red
- ☐ Manatee, Red

### Size

- 0
- 1
- 2
- 3
- 1
- 2
- 3

### Part No.

- |        |
|--------|
| 882005 |
| 882015 |
| 882025 |
| 882035 |
| 882014 |
| 882024 |
| 882034 |

### Foot & Leg Positioning



- ☐ Pommel, Blue
- ☐ Pommel, Blue
- ☐ Pommel, Red
- ☐ Pommel, Red
- ☐ Leg Straps, Blue
- ☐ Leg Straps, Blue
- ☐ Leg Straps, Red
- ☐ Leg Straps, Red

### Size

- 0, 1, 2
- 3
- 1, 2
- 3
- 0, 1, 2
- 3
- 1, 2
- 3

### Part No.

- |        |
|--------|
| 883045 |
| 883095 |
| 883044 |
| 883094 |
| 883055 |
| 883105 |
| 883054 |
| 883104 |

### Trunk Positioning



- ☐ Torso Belt, Blue
- ☐ Torso Belt, Blue
- ☐ Torso Belt, Red
- ☐ Torso Belt, Red

### Size

- 0, 1, 2
- 3
- 1, 2
- 3

### Part No.

- |        |
|--------|
| 883025 |
| 883075 |
| 883024 |
| 883074 |

### Head Support



- ☐ Pads for Head/Lateral Support, Blue
- ☐ Pads for Head/Lateral Support, Blue
- ☐ Pads for Head/Lateral Support, Red
- ☐ Pads for Head/Lateral Support, Red

### Size

- 0, 1, 2
- 3
- 1, 2
- 3

### Part No.

- |        |
|--------|
| 883015 |
| 883065 |
| 883014 |
| 883064 |



### Frame Options

- ☐ Rolling Shower
- ☐ Wall Mount for Frame
- ☐ Frame for Manatee - Tub Stand - 10.5"
- ☐ Tub Slider without Manatee
- ☐ Tub Slider without Manatee
- ☐ Tub Slider with Manatee, Blue (W: 50")
- ☐ Tub Slider with Manatee, Blue (W: 50")
- ☐ Tub Slider with Manatee, Blue (W: 50")
- ☐ Tub Slider with Manatee, Red (W: 50")
- ☐ Tub Slider with Manatee, Red (W: 50")
- ☐ Tub Slider with Manatee, Blue (W: 30")
- ☐ Tub Slider with Manatee, Blue (W: 30")\*
- ☐ Tub Slider with Manatee, Blue (W: 30")\*
- ☐ Tub Slider with Manatee, Red (W: 30")\*
- ☐ Tub Slider with Manatee, Red (W: 30")\*

\* Height adjustable legs, no casters

### Size

- One Size
- One Size
- One Size
- 1, 2
- 3
- 0
- 1
- 2
- 1
- 2
- 0
- 1
- 2
- 1
- 2

### Part No.

882100
882110
882101
MC6000PED
MC6200PED
MC6000PED-1
MC6000PED-2B
MC6000PED-3B
MC6000PED-2R
MC6000PED-3R
5100TXPED-1
5100TXPED-2B
5100TXPED-3B
5100TXPED-2R
5100TXPED-3R

### Accessories

- ☐ Storage Net for Frame, Blue
- ☐ Storage Net for Frame, Red

### Size

- One Size
- One Size

### Part No.

882105
882104

**For shipping rates, please call 1-844-876-6245**


PO Box 1739 Matthews, NC 28106 • 1-844-876-6245 • Fax: 704-882-0751

Effective: 3/22/2017

www.r82.com

Rev 1

## Configuration



A (Back height)  
 B (Arm height)  
 C (Lower leg length)  
 D (Seat depth)  
 E (Chest width)  
 F (Seat width)  
 G (Height from seat to top of head)  
 H (Shoulder width)  
 Age  
 Weight

---

---

---

---

---

---


---

---

---

---

## Measurements



### Manatee

	Size 0 Inches	Size 1 Inches	Size 2 Inches
Seat width (A)	12¾"	12¾"	12¾"
Seat depth (B)	8½ - 9¾"	11¼ - 14½"	15¼ - 18"
Back height (C)	16¾"	16¾"	22½"
Lower leg length (D)	11"	11"	13¼"
Seat height (E)	5½ - 16¾"	5½ - 16¾"	5½ - 16¾"
Total width (G)	16½"	16½"	16½"
Length folded	"40 - 43"	43 - 45¾"	53 - 55¾"
Height folded	5"	5"	5"
Weight	12 lbs.	12 lbs.	13.25 lbs.
Max. load	110 lbs.	110 lbs.	165 lbs.

### Frames

	Shower Frame	Bath Tub Frame
Width	23"	17"
Length	38¼"	28½ - 43¼"
Height	-	10½"
Height, minimum	14¾"	-
Height, minimum, with chair	29¾"	-
Height, middle	17"	-
Height, middle, with chair	32¾"	-
Height, maximum	18¾"	-
Height, maximum, with chair	35¼"	-
Height, folded	8"	-
Weight	17.5 lbs.	4.5 lbs.
Max. load	165 lbs.	165 lbs.