



(Shown with Gunstock Walnut Top Finish and U Base)



Non-Tilt Pneumatic Cylinder Overbed Table

- Pneumatic cylinder lift mechanism – infinite height adjustment from 27" - 42".
- ¾" x 18" x 32" high pressure laminate top with T-molding.
- 2" casters for easy maneuvering.
- Column can be mounted on left or right side.
- Available with "H" (Figure A) or "U" base in chrome or opal powder-coat.
- Optional dual-mirror vanity drawer.
- Base and columns available in additional powder coat finishes.
- Available on drop ship basis only.



A

16008 H

TOP FINISH		H Base Chrome Frame		H Base Opal Frame		UOM
		W/O VANITY	W/VANITY	W/O VANITY	W/VANITY	
Gunstock Walnut		16008-HCGW	16008-HCGWV	16008-HOGW	16008-HOGWV	1/cs
Fusion Maple		16008-HCFM	16008-HCFMV	16008-HOFM	16008-HOFMV	1/cs
Banister Oak		16008-HCBO	16008-HCBOV	16008-HOBO	16008-HOBOV	1/cs
Wild Cherry		16008-HCWC	16008-HCWCV	16008-HOWC	16008-HOWCV	1/cs
English Oak		16008-HCEO	16008-HCEOV	16008-HOEO	16008-HOEOV	1/cs
Figured Mahogany		16008-HCDM	16008-HCDMV	16008-HODM	16008-HODMV	1/cs
Nepal Teak		16008-HCNT	16008-HCNTV	16008-HONT	16008-HONTV	1/cs
Quartered Oak		16008-HCQO	16008-HCQOV	16008-HOQO	16008-HOQOV	1/cs
Solar Oak		16008-HCSO	16008-HCSOV	16008-HOSO	16008-HOSOV	1/cs

16008 U

TOP FINISH		U Base Chrome Frame		U Base Opal Frame		UOM
		W/O VANITY	W/VANITY	W/O VANITY	W/VANITY	
Gunstock Walnut		16008-UCGW	16008-UCGWV	16008-UOGW	16008-UOGWV	1/cs
Fusion Maple		16008-UCFM	16008-UCFMV	16008-UOFM	16008-UOFMV	1/cs
Banister Oak		16008-UCBO	16008-UCBOV	16008-UOBO	16008-UOBOV	1/cs
Wild Cherry		16008-UCWC	16008-UCWCV	16008-UOWC	16008-UOWCV	1/cs
English Oak		16008-UCOE	16008-UCEOV	16008-UOEO	16008-UOEOV	1/cs
Figured Mahogany		16008-UCDM	16008-UCDMV	16008-UODM	16008-UODMV	1/cs
Nepal Teak		16008-UCNT	16008-UCNTV	16008-UONT	16008-UONTV	1/cs
Quartered Oak		16008-UCQO	16008-UCQOV	16008-UOQO	16008-UOQOV	1/cs
Solar Oak		16008-UCSO	16008-UCSOV	16008-UOSO	16008-UOSOV	1/cs

specifications

Height (adjustable)	27" - 42"
Weight	55 lbs.
Carton Shipping Weight	59 lbs.
Weight Limit	100 lbs.
Base Dimensions	18" x 32"
Tabletop Dimensions	18" x 32"
Warranty	Limited Lifetime