



Model no. 12230-35-01



Model no. 12230-02-35-87



Model no. 12220-35-01

Dolomite Alpha Advanced/Basic Rollators

Basket and food & beverage tray standard.

Model no. 12220-35-01 ALPHA Advanced with forearm attachment
Model no. 12230-35-01 ALPHA Basic with forearm attachment

Model no. 12230-02-35-87 ALPHA Basic with platform attachment
Model no. 12220-02-35-87 ALPHA Advanced with platform attachment

Colour: Navy Blue

The Dolomite Alpha is a special rollator designed specifically for users who deal with the added complexity of walking with arthritic or stroke conditions. The unique, adjustable forearm supports or platform option allow for maximum stability. Alpha also offers two braking systems to accommodate the specific needs of the user.

Features

- Frame construction maximizes stability and encourages proper posture and balance
- Forearm supports adjustments in angle, width and depth
- Individually adaptable hand grip can be adjusted in height and angle
- Unique advanced brake grip is extremely easy to use and accommodates limited hand dexterity
- Flip-up seat for greater reach access and gait training
- Dolomite Advanced Brake System (DABS II)

Specifications

	Forearm/Platform
Distance between rear wheels	22"
Minimum handle height	37 ½"/38 ½"
Maximum handle height	47 ½"/48 ½"
Seat height	24 ½"
Overall width	26"
Overall length	27 ½"
Width between handles	14 - 17"
Rollator weight (lbs.)	28 ½"/29 ½"
Weight capacity (lbs.)	275