


Around water, it's safety first

Aeromesh® Rapid-Dry Bath Belt

A secure, comfortable belt is an essential safety feature for any bath or shower chair. Our belt is soft against bare chest, abdomen, knees and legs

- Allows a complete wash and quick dry as water, soap and air flow freely through Aeromesh® fabric
- Quickly open, close and adjust oversized side release buckles and D-Rings



Water and air easily pass through Aeromesh® fabric



Two-piece pad, grommets

Grommet Mount: Screw corrosion-resistant stainless steel grommets to bath chair frame with a screw, Bath Belt Release Knob (shown, right) or Bodypoint mounting clamps (see page 7)



Quick release knobs facilitate cleaning and adjustments (FA812-BP-2, sold separately, in pairs)



One-piece pad, grommets



One-piece pad, plastic slides

Wrap webbing end around the bath chair frame and secure with plastic slides



Universal Elastic Strap

Temporary positioning for work, play, therapy or fittings with highly stretchable, strong neoprene material

- Easy-grasp hook-and-loop end tab with a convenient thumb hole for users with limited hand function
- Use on the chest, abdomen, lap or legs
- Length can be cut to fit

SP110S 142cm (56") x 8cm (3")
SP110L 142cm (56") x 11cm (4-1/2")
SP110XL 142cm (56") x 15cm (6")



Elastic Extension Strap

Combine with the Universal Elastic Strap to increase its length – a solution to fit around large chairs or car seats.

SP111S 72.5cm (28-1/2") x 8cm (3")
SP111L 72.5cm (28-1/2") x 11cm (4-1/2")
SP111XL 72.5cm (28-1/2") x 15cm (6")



ANDREAS

Andreas holds, lifts and plays with Ellie, his new daughter, using the comfort and support of our Universal Elastic Strap.

Learn more about Andreas at bodypoint.com/success-stories



Leg Harness

The "go-to solution" for users that tend to slide out of their chairs, thrust or have spasticity. The padded strap securely holds the leg at the inner thigh without limiting anterior movement or pressing on the bladder. Can be used with or without a hip belt.

LH600S-B1 OR -B2

23cm (9") x 4cm (1-1/2")

LH600M-B1 OR -B2

33cm (13") x 4.5cm (1-3/4")

LH600L-B1 OR -B2

38cm (15") x 5cm (2")



Quick Release Connector

Attach this bracket to any rigid accessory on the wheelchair frame (such as lateral support, headrest mounting bar, or computer support bracket). It locks securely in place without wiggle, yet allows instant release with one hand. Its slim 1/2" profile provides close-fitting attachment along the seat and backrest.

HW311-2