



Human Care’s Convertible Chair is a no-lifting multi-functional device that helps caregivers to more efficiently and effectively perform their safe patient handling practices, such as transfers, repositioning, transport and early mobilization of patients that may help to reduce the patient length of stay, improve patient satisfaction, safety and dignity, and reduce work-related injuries.

I-400

Convertible chair for early mobilization, safe and comfortable transfers

- Unlimited repositioning capabilities
- Central brakes and steer locking functions
- Ideal for acute care setting

The I-400 chair is a safe and comfortable device ideal for acute care settings, although it can just as easily function in a home care or nursing facility setting. It is a key element in early patient mobilization with the added benefits of caregiver safety. Its Tilt-in-Space capability is ideal for the various complications, such as pressure ulcers, experienced by patients with compromised mobility.

The I-400 enables lateral transfer in a supine position safely and easily which builds confidence in both the patient and the caregiver. Its portability allows usage in multiple rooms in a department and the ability to utilize the Patient Transfer System (PTS) simplifies the transfer of patients from a bed to a chair in addition to helping protect the caregiver from injury.



Central locking/steering casters safe for patient and caregiver



Can be equipped with a patient transfer system. See page 10.



Description	Part no.	Seat width	Overall width	Overall length	Overall height	Maximum user weight
I-400 Chair	B1041-TS-GEL	53cm / 21"	76cm / 30"	Prone: 184-199cm / 72.5-78.5" Wheelchair: 115-159cm / 45.5-62.5"	Prone: 70cm / 27.5" Wheelchair: 130-142cm / 51-56"	181 kg / 400 lbs