

EVOLUTION SERIES.

The Evolution Series, available in four sizes, features the strongest and one of the most stable walkers on the market today.

Every Evolution model uses 1-1/8" cold rolled steel tubing, the same material used in mountain bike fabrication. Strong 1" tubing is used on the handles, every fastener is stainless steel or plated against corrosion. The Evolution series offers unmatched stability for users suffering from neurological conditions such as Multiple Sclerosis and Parkinson's, in addition to arthritic and stroke conditions as well as users with wide gate.

FEATURES.

UV Protection Frame Scratch resistant paint, stable footprint

Seat Solid plastic base with removable fitted waterproof cushion

Wire Basket Large & removable




Wheels 8" 40 mm wide, Solid grip / Longer lasting polyurethane, Heavy duty

Plastic Components Withstanding up to -40C

Colour Purplish Red, Sapphire Blue, Silver, Black



SPECIFICATIONS.

LOW	REGULAR	LOW WIDE	WIDE
			
20 lbs 17" Seat Width 19" Seat Height 21" Width x 30" Length 30" - 36" Handle Height 350 lbs Weight Capacity	23 lbs 17" Seat Width 23" Seat Height 21" Width x 30" Length 30" - 38" Handle Height 350 lbs Weight Capacity	22 lbs 21" Seat Width 21" Seat Height 26" Width x 30" Length 30" - 36" Handle Height 400 lbs Weight Capacity	23 lbs 21" Seat Width 23" Seat Height 26" Width x 30" Length 30" - 38" Handle Height 400 lbs Weight Capacity

OPTIONAL ACCESSORIES.

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes
- Flexible Back Strap
- Oxygen Tank Holder
- Mono Brakes
- Tote Bag (fits inside basket)
- Tray - Seat Mounted
- Platform Attachment



Evolution
TECHNOLOGIES*

EVOLUTION TECHNOLOGIES INC.
2530 Davies Avenue,
Port Coquitlam, BC V3C2J9

tf +1 800.556.2558
info@evolution.ca
www.evolution.ca

EVOLUTION TECHNOLOGIES USA INC.
2657 Delta Ring Road
Ferndale WA 98248

tf +1 800.839.2558
info@evolutionwalker.com
www.evolutionwalker.com