



Base Model #201 shown  
Wide Base for Shower Stall use  
Overall Footprint – tallest setting, 20 ½" x 20 ½"

### Specifications 201/202

Adjustable legs vary seat height 19 1/2" – 22 1/2"  
Seat width 18"  
Seat depth 18"  
Weight Capacity 350 lbs



Base Model #202 shown  
Narrow Base for Tub use  
Overall Footprint – crutch tips, 18" x 18"

The name hints at what these products are intended to do. They can be an adjustable height toilet seat, a toilet safety frame, a shower bench, bathtub bench or a transfer bench. The modular design lets you decide.

The densely padded seat has a cut-out that can be positioned to the front, back, or either side and the foam is configured in our unique "**Waterfall**" style that addresses pressure problems caused from sensitive tissue coming in contact with unforgiving surfaces. The open frame system not only affords unimpeded directional seat placement but also allows the optional back and arms to be correctly positioned without restriction and **without tools**. The 350-pound capacity frame is created from high-quality **Stainless Steel** that provides a strong, durable, non-corrosive product.

- ☐ **Model 201 - wide footprint**
- ☐ **Model 202 - narrow footprint**
- ☐ Option A2 One Padded Arm
- ☐ Option A3 Two Padded Arms
- ☐ Option C2 Flat Comfortuff™ Seat
- ☐ Option D7 Padded Back
- ☐ Option K Suction Cups
- ☐ Option BO B-O Transfer Bench
- ☐ Option B Short Transfer Bench
- ☐ Option P Pan Hanger