

# NON-POWERED

Mattress Replacement Systems



## FUSION 1K™

Mattress Replacement System  
Featuring Dynamic Dispersion®

**Manufactured by MOXI Enterprises LLC**

**The first non-powered therapy mattress that converts to true low air loss therapy.**

Dynamic Dispersion differs from all other non-powered alternatives. Each air and foam-filled chamber communicates to create pressure redistribution.

**Size:** FUS1K80: 80" x 36" x 6" (203.2 cm x 91.4 cm x 15.2 cm)  
FUS1K84: 84" x 36" x 6" (213.4 cm x 91.4 cm x 15.2 cm)

**Approximate Weight:** Approximately 46 lbs. (20.9 kg)

**Weight Limit:** 500 lbs. (227 kg) for system

**Warranty:** 5-year non-prorated warranty;  
1-year warranty on cover



**Medicare Code:**  
E0373 Group 2 Support Surfaces

Distributed by ROHO in the USA only

## PRODUCT FEATURES

- The “closed” system envelops and redistributes pressure similar to other cellular-based technologies.
- High resilient foam is encapsulated in preset air chambers with airflow control. Individual cells communicate to evenly distribute the patient’s pressure over the therapy zone.
- Since air, not foam, supports the majority of the load, the system is designed to last years beyond the competition.
- Easily convert your FUSION 1K into a FUSION 2K for full “hospital grade” true low air loss therapy.
- The top stretch cover is anti-microbial, fluid-resistant and breathable.



## FUSION XC™

Dynamic Dispersion®  
Mattress Replacement System

**Manufactured by MOXI Enterprises LLC**

**A simple solution for non-powered therapy.**

Soft, sloped multi-density foam cradles heels to reduce pressure and shearing in critical zones. Each air- and foam-filled chamber communicates to create pressure redistribution.

**Size:** FUSXC80: 80" x 36" x 6" (203.2 cm x 91.4 cm x 15.2 cm)  
FUSXC84: 84" x 36" x 6" (213.4 cm x 91.4 cm x 15.2 cm)

**System weight:** Approximately 46 lbs. (20.9 kg)

**Weight limit:** 500 lbs. (227 kg)

**Warranty:** 5-year limited warranty on the mattress; 1-year limited warranty on the cover



**Medicare Code:**  
E0373 Group 2 Support Surfaces

Distributed by ROHO in the USA only

## PRODUCT FEATURES

- The system, immerses and envelops in order to redistribute pressure.
- Top stretch cover is anti-microbial, fluid-resistant and breathable.
- Individual chambers communicate to evenly distribute the patient’s pressure over the therapy zone.
- Since air, not foam, supports the majority of the load, the system is designed to last years beyond the competition.
- High resilient foam is encapsulated in preset air chambers with airflow control. Individual cells communicate to evenly distribute the patient’s pressure over the therapy zone.