

## ROHO® MINI-MAX® Cushion

This is a cushion designed for short-term use during recreational and sport activities, bathing, traveling, etc. The MINI-MAX provides shock and impact absorption, as well as skin protection. The 1.25 in. (3.2 cm) air-filled cells are fully adjustable.

- ✓ Comfortable
- ✓ Absorbs shock
- ✓ Adjustable

### Dimensions

Height: 1.25 in. / 3 cm  
 Width: 15 in. to 18.25 in. /  
 38.5 cm to 46.5 cm  
 Depth: 14.75 in. to 16 in. /  
 37.5 cm to 40.5 cm  
 Weight: 2 lbs. / 0.9 kg

### Weight limit

None, if properly sized

### Limited warranty

Cushion: 24 months  
 Cover: 6 months



## ROHO® PACK-IT® Cushion

The ROHO PACK-IT Cushion is a smaller version of the MINI-MAX. It's designed for athletic competitions, bathing, traveling or as a scapular, sacral or lumbar support. The PACK-IT's size and built-in strap makes it easy to attach and remove the cushion on the go.

- ✓ Strap for easy attachment and removal
- ✓ Absorbs shock
- ✓ Skin/soft tissue protection
- ✓ Adjustable

### Dimensions

Height: 1.25 in. / 3 cm  
 Width: 16.25 in. / 41.5 cm  
 Depth: 9.5 in. / 24 cm  
 Weight: 1.5 lbs. / 0.7 kg

### Weight limit

None, if properly sized.

### Limited warranty

Cushion: 24 months  
 Cover: 6 months

